

ALUMA ELITE MAT TABLE

STANDARD FEATURES:

- Load capacity: 1000 lbs
- Durable aluminum construction
- Tough powder coat finish
- Leveling feet
- 1.5" Semi-firm foam cushion
- Rounded corners
- Elevation option: hand control is mounted on one side
- Graphics package: includes (1) school logo/name decal on each tape holder (1 on each end) (Only available with storage shelf option)



A4520 with
Storage Shelves

STANDARD FIXED HEIGHT MAT TABLES:

A4495

48"W x 84"L x 21"H (4 Legs)



A4520

60"W x 84"L x 21"H (4 Legs)



A4525

72"W x 96"L x 21"H (4 Legs)



STANDARD ELEVATING MAT TABLES:

A4495-EL

48"W x 84"L x 21-31"H (4 Legs)




A4520-EL

60"W x 84"L x 21-31"H (4 Legs)



A4525-EL

72"W x 96"L x 21-31"H (4 Legs)



OPTIONS:



A4520-EL with Storage Shelves

OPTIONS:

A4495A - Storage Shelf for A4495 & A4495-EL (One on Each End)

A4520A - Storage Shelf for A4520 & A4520-EL (One on Each End)

A4525A - Storage Shelf for A4525 & A4525-EL (One on Each End)

S1075AE - Single Foot Control Pedal (Elevating Mat Table Only)

SBO-047 - Side Bumper for Mat Table, 4'W x 7'L (Not available with A4495A, A4520A, or A4525A)

SBO-057 - Side Bumper for Mat Table, 5'W x 7'L (Not available with A4495A, A4520A, or A4525A)

SBO-068 - Side Bumper for Mat Table, 6'W x 8'L (Not available with A4495A, A4520A, or A4525A)

STANDARD FRAME COLORS:

	Black		Burgundy
	Light Gray		Green
	Charcoal		Forest Green
	Bright Blue		Purple
	Royal Blue		Yellow
	Navy Blue		Orange
	Red		Brown
	White		Sky Blue

ACCENT COLOR FOR STORAGE SHELF:

	Black		Burgundy
	Light Gray		Green
	Charcoal		Forest Green
	Bright Blue		Purple
	Royal Blue		Yellow
	Navy Blue		Orange
	Red		Brown
	White		Sky Blue

STANDARD UPHOLSTERY CHOICES:

	Royal Blue		American Beauty
	Imperial Blue		Lark
	Blue Ridge		Burgundy
	Deep Violet		Dove
	Forest		Charcoal
	Mandarin Orange		Mocha
	Black		Putty