

PT2000 BARIATRIC HI-LO MAT TABLE

STANDARD FEATURES:

- Lifting capacity: 1000 lbs.
- Height range: 18-33"
- Durable steel base with powder coat finish
- 1.5" High density foam cushion
- Rounded corners
- Single foot control
- Convenient drop down casters
- ADA Compliant



PT2000-47



PT2000-47

STANDARD TABLE:

PT2000-47

48"W x 84"L x 18-33"H

PT2000-57

60"W x 84"L x 18-33"H

PT2000-68

72"W x 96"L x 18-33"H

OPTIONS:

SBO-047 - Side Bumper for Mat Table, 4'W x 7'L

SBO-057 - Side Bumper for Mat Table, 5'W x 7'L

SBO-068 - Side Bumper for Mat Table, 6'W x 8'L



STANDARD UPHOLSTERY CHOICES:

	Royal Blue		American Beauty
	Imperial Blue		Lark
	Blue Ridge		Burgundy
	Deep Violet		Dove
	Forest		Charcoal
	Mandarin Orange		Mocha
	Black		Putty