

PT2001 BARIATRIC HI-LO MAT TABLE

STANDARD FEATURES:

- Lifting capacity: 750 lbs.
- Height range: 18-33"
- Durable steel base with powder coat finish
- 1.5" High density foam cushion
- Rounded corners
- Single foot control
- ADA Compliant



PT2001-47

STANDARD TABLE:

PT2001-47

48"W x 84"L x 18-33"H

PT2001-57

60"W x 84"L x 18-33"H

PT2001-68

72"W x 96"L x 18-33"H

OPTIONS:

SBO-047 - Side Bumper for Mat Table, 4"W x 7"L

SBO-057 - Side Bumper for Mat Table, 5"W x 7"L

SBO-068 - Side Bumper for Mat Table, 6"W x 8"L



STANDARD UPHOLSTERY CHOICES:

	Royal Blue		American Beauty
	Imperial Blue		Lark
	Blue Ridge		Burgundy
	Deep Violet		Dove
	Forest		Charcoal
	Mandarin Orange		Mocha
	Black		Putty