

■ PREMIER COLLECTION

# Ascent

Available: Power, Manual, Power XL, Manual XL

The Ascent is an upscale, stylish recliner, featuring swing-away arms for patient transport and extra padding for increased comfort. Also available in wingback style, the Ascent boasts the durable engineering and strength of the Classic but with an upscale design, making it an ideal medical chair for patients and providers.



## Key Benefits

- With a reinforced steel powder-coated frame and heavy-duty casters, the Ascent medical recliner is built to last.
- Designed to accommodate patients up to 500 pounds, the bariatric model (XL, 66P) is deeper and wider than the standard model.
- For quality engineering and patient-centric innovation that prioritizes functionality, safety, and comfort, Champion products comply with ISO 13485 and are made in the USA.

## Common Applications:

Oncology & Infusion Settings

## Standard Features



Trendelenburg Positioning



Foot Tray



Fold-Down Side Tablet



4" (10.2 cm) Casters




Push Handle


(Additional standard features are also included. Talk to a Sales Rep for details.)

## Space Requirements


Front View




Top View



Side View



Top View with Swing Away Arms



Shown with optional second, fold away table.


## Specifications

	Standard Model	XL Model
<b>Power Recline</b>	Yes, Model: 65P	Yes, Model: 66P
<b>Manual Recline</b>	Yes, Model: 650	Yes, Model: 660
<b>Weight Capacity</b>	300 lbs. / 136 kg	500 lbs. / 226 kg
<b>Back Height</b> (from floor ♦ from seat)	47" ♦ 29" 119.40 cm ♦ 73.70 cm	48" ♦ 30" 121.92 cm ♦ 76.20 cm
<b>Overall Width</b> (with 2 gray tables down ♦ up)	32.25" ♦ 52" 81.92 cm ♦ 132.10 cm	36" ♦ 56" 91.44 cm ♦ 142.24 cm
<b>Maximum Width</b> (arms open at widest)	90" / 228.60 cm	94" / 238.76 cm
<b>Minimum Width</b> (arms open at 180°)	32.25" / 81.920 cm	36.50" / 92.71 cm
<b>Seat Width</b>	22.75" / 57.79 cm	26.75" / 67.95 cm
<b>Seat Height ♦ Seat Depth</b>	24.50" ♦ 20.50" 62.23 cm ♦ 52.32 cm	24.50" ♦ 20.50" 62.23 cm ♦ 52.32 cm
<b>Arm Height</b> (highest point from floor)	29" / 73.70 cm	29" / 73.70 cm
<b>Overall Upright Depth ♦ Trendelenburg Depth</b>	43" ♦ 75" 109.22 cm ♦ 190.50 cm	43" ♦ 76" 109.22 cm ♦ 193.04 cm
<b>Dual Thermoplastic Casters</b> (diameter)	4" / 10.20 cm	4" / 10.20 cm

## Available Options



Heat + Massage Tabletop Switch



USB Charging Port



Fold-Away Table With Cupholder



Paper Roll Dispenser



Central-Locking Caster System