



EDGE SPORT WOOD MAT TABLE



STANDARD FEATURES:

- Load capacity: 450 lbs.
- Solid hardwood frame and 4 legs
- 3/4" Plywood top
- 1.5" High density foam cushion
- Rounded corners
- Triple-bolted leg assembly
- ADA Compliant

STANDARD TABLES:

MAT-4872-4L 48"W x 72"L x 19"H	MAT-4884-4L 48"W x 84"L x 19"H	MAT-6084-4L 60"W x 84"L x 19"H

OPTIONS:










LOGO-BB - Table Sling





EDGE SPORT WOOD MAT TABLE

STANDARD UPHOLSTERY COLORS:

 Royal Blue	 American Beauty
 Imperial Blue	 Lark
 Blue Ridge	 Burgundy
 Deep Violet	 Dove
 Forest	 Charcoal
 Mandarin Orange	 Mocha
 Black	 Putty








STANDARD WOOD COLORS:

 Oak
 Maple

WOOD STAIN OPTIONS:

MAT-Stain

STANDARD STAIN COLORS:

 Warm Copper	 Rustic Fall
 Deep Garnet	 Warm Caramel
 Fossil Gray	 Black Beauty
 Tungsten	Custom Wood Stain Color: