



# 6 LEG EDGE

# SPORT WOOD MAT TABLE



MAT-4884-6L

### STANDARD FEATURES:

- Load capacity: 700 lbs.
- Solid hardwood frame and 6 legs
- 3/4" Plywood top
- 1.5" High density foam cushion
- Rounded corners
- Triple-bolted leg assembly
- ADA Compliant

### STANDARD TABLES:

<input type="radio"/> MAT-4884-6L 48"W x 84"L x 19"H	<input type="radio"/> MAT-6084-6L 60"W x 84"L x 19"H	<input type="radio"/> MAT-7296-6L 72"W x 96"L x 19"H
QTY:	QTY:	QTY:

### OPTIONS:

- LOGO-BB - Table Sling





# 6 LEG EDGE SPORT WOOD MAT TABLE

## STANDARD UPHOLSTERY COLORS:

- |   |   |
|---|---|
| <input type="radio"/>  Royal Blue      | <input type="radio"/>  American Beauty |
| <input type="radio"/>  Imperial Blue   | <input type="radio"/>  Lark            |
| <input type="radio"/>  Blue Ridge      | <input type="radio"/>  Burgundy        |
| <input type="radio"/>  Deep Violet     | <input type="radio"/>  Dove            |
| <input type="radio"/>  Forest          | <input type="radio"/>  Charcoal        |
| <input type="radio"/>  Mandarin Orange | <input type="radio"/>  Mocha           |
| <input type="radio"/>  Black           | <input type="radio"/>  Putty           |







## STANDARD WOOD COLORS:

- |   |
|---|
| <input type="radio"/>  Oak    |
| <input type="radio"/>  Maple |

## WOOD STAIN OPTIONS:

- |                                 |
|---------------------------------|
| <input type="radio"/> MAT-Stain |
|---------------------------------|

## STANDARD STAIN COLORS:

- |   |  |
|---|--|
| <input type="radio"/>  Warm Copper | <input type="radio"/>  Rustic Fall  |
| <input type="radio"/>  Deep Garnet | <input type="radio"/>  Warm Caramel |
| <input type="radio"/>  Fossil Gray | <input type="radio"/>  Black Beauty |
| <input type="radio"/>  Tungsten    | <input type="radio"/> Custom Wood Stain Color:   |