



The Wall Mount Mat Table is a fold out mat table ideal for small spaces. The table is constructed of solid hardwood with a flat 1.5" high density foam cushion with softly rounded corners. A platform attaches to the wall in the back, while the three front legs and cushion fold out away from the wall when needed. The triple-bolted leg assembly adds extra durability.

STANDARD FEATURES:

- Load Capacity: 450 lbs.
- Solid hardwood frame and legs
- 3/4" Plywood top
- 1.5" High density foam flat cushion
- Rounded corners
- Triple-bolted leg assembly
- 10" Depth when folded
- ADA Compliant

Wall Mount
MAT TABLES




WMMT-4884

STANDARD TABLES:

WMMT-4884
48"W x 84"L x 19"H

WMMT-6084
60"W x 84"L x 19"H



STANDARD WOOD CHOICES:



Oak



Maple

WOOD STAIN OPTION:

MAT-Stain

STAIN CHOICES:

- Warm Copper
- Deep Garnet
- Fossil Gray
- Tungsten

- Rustic Fall
- Warm Caramel
- Black Beauty
- Custom Wood Stain Color

STANDARD UPHOLSTERY CHOICES:

- | | |
|-----------------|-----------------|
| Royal Blue | American Beauty |
| Imperial Blue | Lark |
| Blue Ridge | Burgundy |
| Deep Violet | Dove |
| Forest | Charcoal |
| Black | Mocha |
| Mandarin Orange | Putty |