

■ PREMIER COLLECTION

Ascent II

Available: Power, Power XL

Keeping patients comfortable and safe during long infusion or oncology treatments is crucial. With the Ascent II, we ensure that healthcare professionals have clinical access and ease of use, while creating a comforting and reassuring experience for the patient.



Key Benefits

- As an extension of the Ascent line, the Ascent II recliner features a “wingback” style that offers additional headroom for the patient and unique two-tone customization options.
- Swing-away arms and a removable seat aid in meeting infection control goals.
- The four-inch dual thermoplastic casters are corrosion resistant and roll with ease.
- For quality engineering and patient-centric innovation that prioritizes functionality, safety, and comfort, Champion products comply with ISO 13485 and are made in the USA.

Common Applications:

Oncology & Infusion Settings

Standard Features



Trendelenburg Positioning



Swing-Away Arms



Foot Tray



Fold-Down Side Tablet




Push Handle


(Additional standard features are also included. Talk to a Sales Rep for details.)

Space Requirements


Front View




Top View



Side View



Top View with Swing Away Arms



Shown with optional second, fold away table.


Specifications

	POWER 67P	Power XL 68P
Power Recline	Yes, Model: 67P	Yes, Model: 68P
Weight Capacity	300 lbs. / 136 kg	500 lbs. / 226 kg
Back Height (from floor ♦ from seat)	47" ♦ 29" 119.40 cm ♦ 73.70 cm	48" ♦ 30" 121.92 cm ♦ 76.20 cm
Overall Width (with 2 gray tables down ♦ up)	32.25" ♦ 52" 81.92 cm ♦ 132.10 cm	36" ♦ 56" 91.44 cm ♦ 142.24 cm
Maximum Width (arms open at widest)	90" / 228.60 cm	94" / 238.76 cm
Minimum Width (arms open at 180°)	32.25" / 81.92 cm	36.50" / 92.71 cm
Seat Width	22.75" / 57.79 cm	26.75" / 67.95 cm
Seat Height ♦ Seat Depth	24.50" ♦ 20.50" 62.23 cm ♦ 52.32 cm	24.50" ♦ 20.50" 62.23 cm ♦ 52.32 cm
Arm Height (highest point from floor)	29" / 73.70 cm	29" / 73.70 cm
Overall Upright Depth ♦ Trendelenburg Depth	43" ♦ 75" 109.22 cm ♦ 190.50 cm	43" ♦ 76" 109.22 cm ♦ 193.04 cm
Dual Thermoplastic Casters (diameter)	4" / 10.20 cm	4" / 10.20 cm


Available Options



Stylish Black Arm Caps



USB Charging Port



Paper Roll Dispenser



Central-Locking Caster System



Heat + Massage Tabletop Switch