


# VT65MLGP



33.5" x 23.63" x 23.5" (H x W x D)

Freestanding medical general purpose all-freezer capable of -25° C operation with front-mounted lock

## Highlights:

Ideal 24" footprint for easy placement in healthcare, institutional, and similar settings

-25°C operation for safe storage of temperature sensitive materials

Digital thermostat for easy monitoring of interior temperature

## Product Features:

Healthcare freezer	Designed for general purpose storage in institutional facilities
-25°C operation	Low temperature operation is ideal for certain vaccines, pharmaceuticals, and other temperature sensitive materials
Factory installed lock	Keyed lock for a secure interior
Three removable storage drawers	Wire drawers with transparent plastic fronts allow for easy storage organization
Manual defrost	Reduced temperature fluctuations with an energy efficient manual defrost system
Digital thermostat	Digital thermostat for easy monitoring of interior temperature
One-piece interior liner	Easy clean-up with a seamless liner that won't hold a mess in hidden crevices
Reversible door	User reversible door swing for added flexibility
Fully finished white cabinet	Allows the unit to be used freestanding

### VT65MLGP Specifications:

Overview	
Nominal Height	34.0" (86 cm)
Height of Cabinet	33.5" (85 cm)
Nominal Width	24.0" (61 cm)
Width	23.63" (60 cm)
Nominal Depth	24.0" (61 cm)
Depth	23.5" (60 cm)
Capacity	3.5 cu.ft. (99 L)
Defrost Type	Manual
Door	White
Cabinet	White
US Electrical Safety	UL
Voltage/Frequency	115 V AC/60 Hz
Weight	105.0 lbs. (48 kg)
Parts & Labor Warranty	1 Year
Compressor Warranty	5 Years
UPC	761101068152
Freezer	
Door Swing	RHD
Reversible	Yes
Interior Height	24.25" (62 cm)
Interior Width	17.75" (45 cm)
Interior Depth	16.5" (42 cm)
Thermostat Type	Digital
Refrigerant Type	R134a
Refrigerant Amount	2.76 oz.
Level Legs	2
Interior Drawers	Yes
Drawers Qty	3
Temperature Range	-18 to -25°C

