

# VT65MLPLUS2ADA



32.25" x 23.63" x 23.5" (H x W x D)

ADA compliant 24" wide all-freezer for freestanding use, manual defrost with a NIST calibrated thermometer, lock, and -25°C capability

## Highlights:

32" height for use in ADA compliant settings

NIST calibrated thermometer displays the current and high/low temperature to the nearest tenth of a degree

-25° C operation for safe storage of temperature sensitive materials

## Product Features:

ADA compliant	32" height for use in ADA compliant settings
Factory installed lock	Keyed lock for secure interior
NIST calibrated temperature display	Thermometer provides continuous display of the current and high/low temperature to the nearest tenth of a degree, with a readout switchable from Celsius to Fahrenheit at the touch of a button
Buffered temperature probe	Thermometer probe is encased in a glycol-filled bottle to better simulate the temperature of stored contents instead of freezer air temperature
-25°C capable	Low temperature capability ideal for certain sensitive items
Digital thermostat	Electronic controls located in the front kickplate for easy and intuitive temperature management
Manual defrost	Efficient performance with reliable temperature stability
Pull-out drawers	Clear baskets for convenient and separated storage
Magnetic door gasket	Quality door gaskets ensure a tight seal for energy savings and improved performance

## VT65MLPLUS2ADA Specifications:

Overview	
Nominal Height	32.0" (81 cm)
Height of Cabinet	32.0" (81 cm)
Height to Hinge Cap	32.25" (82 cm)
Nominal Width	24.0" (61 cm)
Width	23.63" (60 cm)
Nominal Depth	24.0" (61 cm)
Depth	23.5" (60 cm)
Depth with Handle	24.25" (62 cm)
Depth with door at 90°	45.0" (114 cm)
Capacity	3.2 cu.ft. (91 L)
Defrost Type	Manual
Door	White
Cabinet	White
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Amps	1.3
Voltage/Frequency	115 V AC/60 Hz
Weight	105.0 lbs. (48 kg)
Parts & Labor Warranty	1 Year
Compressor Warranty	5 Years
UPC	761101058054

Freezer	
Door Swing	RHD
Reversible	Factory Reversible
Interior Height	24.25" (62 cm)
Interior Width	17.75" (45 cm)
Interior Depth	16.5" (42 cm)
Thermostat Type	Digital
Refrigerant Type	R134a
Refrigerant Amount	2.7 oz.
High Side PSI	235.0
Low Side PSI	70.0
Interior Drawers	Yes
Drawers Qty	3
Compressor Step Height	6.5" (17 cm)
Compressor Step Width	17.5" (44 cm)
Compressor Step Depth	5.5" (14 cm)

### Temperature Range

