

# Labeling for Cryogenic and Cold Storage

Cold storage labels use a special adhesive that is designed to withstand extreme cold temperatures as low as  $-196^{\circ}\text{C}$  ( $-311^{\circ}\text{F}$ ). Our cold storage labels can be printed or written on and will not bleed, smudge or smear. These labels are offered in a wide variety of sizes and colors on rolls or laser sheets--perfect for distinguishing various specimens after prolonged periods of cold storage.


## Features:

- Able to withstand temperatures from  $-196^{\circ}\text{C}$  up to  $+80^{\circ}\text{C}$
- Can be written on
- Available on rolls or sheets
- Full line of colors
- Various sizes offered
- Label sheets are compatible with laser printers
- Printing instructions and template included




# Cryogenic and Cold Storage Labels


## Label Sheets




LCS-D95W



LCS-D13RW



LCS-24X13W



LCS-33X13W

### Dots for Tubes

Item #	Description	Quantity	Color
LCS-D95W	9.5mm, For .5-1.5ml Tubes	20 Sheets, 192 Labels Per Sheet	White
LCS-D95B	9.5mm, For .5-1.5ml Tubes	20 Sheets, 192 Labels Per Sheet	Blue
LCS-D95G	9.5mm, For .5-1.5ml Tubes	20 Sheets, 192 Labels Per Sheet	Green
LCS-D95N	9.5mm, For .5-1.5ml Tubes	20 Sheets, 192 Labels Per Sheet	Orange
LCS-D95R	9.5mm, For .5-1.5ml Tubes	20 Sheets, 192 Labels Per Sheet	Red
LCS-D95Y	9.5mm, For .5-1.5ml Tubes	20 Sheets, 192 Labels Per Sheet	Yellow
LCS-D95RW	9.5mm, For .5-1.5ml Tubes	20 Sheets, 960 Labels per Color	Assorted
LCS-D13W	13mm, For 1.5-2ml Tubes	20 Sheets, 192 Labels Per Sheet	White
LCS-D13B	13mm, For 1.5-2ml Tubes	20 Sheets, 192 Labels Per Sheet	Blue
LCS-D13G	13mm, For 1.5-2ml Tubes	20 Sheets, 192 Labels Per Sheet	Green
LCS-D13N	13mm, For 1.5-2ml Tubes	20 Sheets, 192 Labels Per Sheet	Orange
LCS-D13R	13mm, For 1.5-2ml Tubes	20 Sheets, 192 Labels Per Sheet	Red
LCS-D13Y	13mm, For 1.5-2ml Tubes	20 Sheets, 192 Labels Per Sheet	Yellow
LCS-D13RW	13mm, For 1.5-2ml Tubes	20 Sheets, 960 Labels Per Color	Assorted

Assorted colors contains Blue, Green, Violet, Red, and Yellow


### Labels for Tubes

Item #	Description	Quantity	Color
LCS-24X13W	24x13mm, For .5 Tubes	20 Sheets, 119 Labels Per Sheet	White
LCS-24X13B	24x13mm, For .5 Tubes	20 Sheets, 119 Labels Per Sheet	Blue
LCS-24X13G	24x13mm, For .5 Tubes	20 Sheets, 119 Labels Per Sheet	Green
LCS-24X13N	24x13mm, For .5 Tubes	20 Sheets, 119 Labels Per Sheet	Orange
LCS-24X13R	24x13mm, For .5 Tubes	20 Sheets, 119 Labels Per Sheet	Red
LCS-24X13Y	24x13mm, For .5 Tubes	20 Sheets, 119 Labels Per Sheet	Yellow
LCS-24X13RW	24x13mm, For .5 Tubes	20 Sheets, 595 Labels Per Color	Assorted
LCS-33X13W	33x13mm, For 1.5-2ml Tubes	20 Sheets, 85 Labels Per Sheet	White
LCS-33X13B	33x13mm, For 1.5-2ml Tubes	20 Sheets, 85 Labels Per Sheet	Blue
LCS-33X13G	33x13mm, For 1.5-2ml Tubes	20 Sheets, 85 Labels Per Sheet	Green
LCS-33X13N	33x13mm, For 1.5-2ml Tubes	20 Sheets, 85 Labels Per Sheet	Orange
LCS-33X13R	33x13mm, For 1.5-2ml Tubes	20 Sheets, 85 Labels Per Sheet	Red
LCS-33X13Y	33x13mm, For 1.5-2ml Tubes	20 Sheets, 85 Labels Per Sheet	Yellow
LCS-33X13RW	33x13mm, For 1.5-2ml Tubes	20 Sheets, 425 Labels per Color	Assorted

Assorted colors contains Blue, Green, Violet, Red, and Yellow

# Cryogenic and Cold Storage Labels

## Label Sheets




LCS-25X25RW

### For Microplates

Item #	Description	Quantity	Color
LCS-38X6W	38x6mm	20 Sheets, 156 Labels Per Sheet	White <input type="checkbox"/>
LCS-38X6B	38x6mm	20 Sheets, 156 Labels Per Sheet	Blue <input type="checkbox"/>
LCS-38X6G	38x6mm	20 Sheets, 156 Labels Per Sheet	Green <input type="checkbox"/>
LCS-38X6N	38x6mm	20 Sheets, 156 Labels Per Sheet	Orange <input type="checkbox"/>
LCS-38X6R	38x6mm	20 Sheets, 156 Labels Per Sheet	Red <input type="checkbox"/>
LCS-38X6Y	38x6mm	20 Sheets, 156 Labels Per Sheet	Yellow <input type="checkbox"/>
LCS-38X6RW	38x6mm	20 Sheets, 780 Labels per Color	Assorted <input type="checkbox"/>
LCS-67X5W	67x5mm	20 Sheets, 99 Labels Per Sheet	White <input type="checkbox"/>


Assorted colors contains Blue, Green, Violet, Red, and Yellow



LCS-25X25W

### For Microscope Slides

Item #	Description	Quantity	Color
LCS-25X25W	25x25mm	20 Sheets, 80 Labels Per Sheet	White <input type="checkbox"/>



LCS-38X19W

### For General Use


Item #	Description	Quantity	Color
LCS-38X19W	38x19mm	20 Sheets, 60 Labels Per Sheet	White <input type="checkbox"/>
LCS-38X19B	38x19mm	20 Sheets, 60 Labels Per Sheet	Blue <input type="checkbox"/>
LCS-38X19G	38x19mm	20 Sheets, 60 Labels Per Sheet	Green <input type="checkbox"/>
LCS-38X19N	38x19mm	20 Sheets, 60 Labels Per Sheet	Orange <input type="checkbox"/>
LCS-38X19R	38x19mm	20 Sheets, 60 Labels Per Sheet	Red <input type="checkbox"/>
LCS-38X19Y	38x19mm	20 Sheets, 60 Labels Per Sheet	Yellow <input type="checkbox"/>
LCS-38X19RW	38x19mm	20 Sheets, 300 Labels per Color	Assorted <input type="checkbox"/>

Assorted colors contains Blue, Green, Violet, Red, and Yellow

# Cryogenic and Cold Storage Labels

## Label Sheets

### For Cryovials




LCS-38X19RW

Item #	Description	Quantity	Color
LCS-43X19W	43x19mm	20 Sheets, 52 Labels Per Sheet,	White <input type="checkbox"/>
LCS-43X19B	43x19mm	20 Sheets, 52 Labels Per Sheet	Blue <input type="checkbox"/>
LCS-43X19G	43x19mm	20 Sheets, 52 Labels Per Sheet	Green <input type="checkbox"/>
LCS-43X19N	43x19mm	20 Sheets, 52 Labels Per Sheet	Orange <input type="checkbox"/>
LCS-43X19R	43x19mm	20 Sheets, 52 Labels Per Sheet	Red <input type="checkbox"/>
LCS-43X19Y	43x19mm	20 Sheets, 52 Labels Per Sheet	Yellow <input type="checkbox"/>
LCS-43X19RW	43x19mm	20 Sheets, 260 Labels Per Color	Assorted <input type="checkbox"/>

Assorted colors contains Blue, Green, Violet, Red, and Yellow

### For Racks and Boxes



LCS-67X25RW

Item #	Description	Quantity	Color
LCS-67X25W	67x25mm	20 Sheets, 30 Labels Per Sheet	White <input type="checkbox"/>
LCS-67X25B	67x25mm	20 Sheets, 30 Labels Per Sheet	Blue <input type="checkbox"/>
LCS-67X25G	67x25mm	20 Sheets, 30 Labels Per Sheet	Green <input type="checkbox"/>
LCS-67X25N	67x25mm	20 Sheets, 30 Labels Per Sheet	Orange <input type="checkbox"/>
LCS-67X25R	67x25mm	20 Sheets, 30 Labels Per Sheet	Red <input type="checkbox"/>
LCS-67X25Y	67x25mm	20 Sheets, 30 Labels Per Sheet	Yellow <input type="checkbox"/>
LCS-67X25RW	67x25mm	20 Sheets, 150 Labels Per Color	Assorted <input type="checkbox"/>
LCS-28X44W	28x44mm	42 Labels Per Sheet	White <input type="checkbox"/>
LCS-64X13W	64x13mm	20 Sheets, 60 Labels Per Sheet	White <input type="checkbox"/>
LCS-76X51W	76x51mm	20 Sheets, 12 labels Per Sheet	White <input type="checkbox"/>

Assorted colors contains Blue, Green, Violet, Red, and Yellow

# Cryogenic and Cold Storage Labels

## Label Rolls



LCR-95W



LCR-95RW




LCR-13R



LCR-19W



LCR-95BIO



LCR-95RAD

### Dots for Tubes

Item #	Description	Quantity	Color
LCR-95W	9.5mm, For .5-1.5ml Tubes	1000 Labels per Roll	White
LCR-95B	9.5mm, For .5-1.5ml Tubes	1000 Labels per Roll	Blue
LCR-95C	9.5mm, For .5-1.5ml Tubes	1000 Labels per Roll	Coffee
LCR-95GL	9.5mm, For .5-1.5ml Tubes	1000 Labels per Roll	Gold
LCR-95G	9.5mm, For .5-1.5ml Tubes	1000 Labels per Roll	Green
LCR-95L	9.5mm, For .5-1.5ml Tubes	1000 Labels per Roll	Lavender
LCR-95N	9.5mm, For .5-1.5ml Tubes	1000 Labels per Roll	Orange
LCR-95P	9.5mm, For .5-1.5ml Tubes	1000 Labels per Roll	Pink
LCR-95R	9.5mm, For .5-1.5ml Tubes	1000 Labels per Roll	Red
LCR-95SL	9.5mm, For .5-1.5ml Tubes	1000 Labels per Roll	Silver
LCR-95Y	9.5mm, For .5-1.5ml Tubes	1000 Labels per Roll	Yellow
LCR-95RW	9.5mm, For .5-1.5ml Tubes	1000 Dots each color	Assorted
LCR-13W	13mm, For 1.5-2ml Tubes	1000 Labels per Roll	White
LCR-13B	13mm, For 1.5-2ml Tubes	1000 Labels per Roll	Blue
LCR-13C	13mm, For 1.5-2ml Tubes	1000 Labels per Roll	Coffee
LCR-13GL	13mm, For 1.5-2ml Tubes	1000 Labels per Roll	Gold
LCR-13G	13mm, For 1.5-2ml Tubes	1000 Labels per Roll	Green
LCR-13L	13mm, For 1.5-2ml Tubes	1000 Labels per Roll	Lavender
LCR-13N	13mm, For 1.5-2ml Tubes	1000 Labels per Roll	Orange
LCR-13P	13mm, For 1.5-2ml Tubes	1000 Labels per Roll	Pink
LCR-13R	13mm, For 1.5-2ml Tubes	1000 Labels per Roll	Red
LCR-13SL	13mm, For 1.5-2ml Tubes	1000 Labels per Roll	Silver
LCR-13Y	13mm, For 1.5-2ml Tubes	1000 Labels per Roll	Yellow
LCR-13RW	13mm, For 1.5-2ml Tubes	1000 Dots each color	Assorted
LCR-19W	19mm Dots, For 15-50ml Tubes	1000 Labels per Roll	White
LCR-21X7W	21x7mm, For 0.5ml Tubes	1500 Labels per Roll	White
LCR-24X13W	24x13mm, For 0.5ml Tubes	1000 Labels per Roll	White
LCR-25W	25mm Dots, For 50ml Tubes	500 Labels per Roll	White

Assorted colors contains Blue, Green, Violet, Red, and Yellow







### Hazard Dots

Item #	Description	Quantity	Color
LCR-95BIO	9.5mm, Biohazard Dots	1000 Labels per Roll	Red/Black
LCR-95RAD	9.5mm, Radioactive Dots	1000 Labels per Roll	Yellow/Pink

# Cryogenic and Cold Storage Labels

## Label Rolls

### For Tubes

Item #	Description	Quantity	Color
LCR-33X13W	33x13mm, For 1.5/2ml Tubes	1000 Labels per Roll	White 
LCR-33X13B	33x13mm, For 1.5/2ml Tubes	1000 Labels per Roll	Blue 
LCR-33X13G	33x13mm, For 1.5/2ml Tubes	1000 Labels per Roll	Green 
LCR-33X13N	33x13mm, For 1.5/2ml Tubes	1000 Labels per Roll	Orange 
LCR-33X13R	33x13mm, For 1.5/2ml Tubes	1000 Labels per Roll	Red 
LCR-33X13Y	33x13mm, For 1.5/2ml Tubes	1000 Labels per Roll	Yellow 

LCR-33X13W



### For General Use

Item #	Description	Quantity	Color
LCR-38X19W	38x19mm	1000 Labels per Roll	White 
LCR-38X19B	38x19mm	1000 Labels per Roll	Blue 
LCR-38X19G	38x19mm	1000 Labels per Roll	Green 
LCR-38X19N	38x19mm	1000 Labels per Roll	Orange 
LCR-38X19R	38x19mm	1000 Labels per Roll	Red 
LCR-38X19Y	38x19mm	1000 Labels per Roll	Yellow 

LCR-38X19R



## Label Dispenser

Item #	Description	Quantity	Color
LCR-DIS	Multi-Roll Dispenser	1 each	Clear

LCR-DIS




Perfectly fits boxes of label rolls for easy dispensing and organization.

# Cryogenic and Cold Storage Labels

## Direct Thermal Label Rolls

**Direct thermal Cryogenic Storage Labels share the same features and benefits as our standard cryo line with one major added perk--they work in your Direct Thermal Printer!**


**With a wide variety of sizes, colors and configurations you are sure to find a solution perfect for your existing setup.**




LTC-95W



LTC-38X13W















LTC-51X6W















LTC-22X22W

### Dots for Tubes

Item #	Description	Quantity	Color
LTC-95W	9.5mm, For .5-1.5ml Tubes	4000 Labels per Roll	White 
LTC-95B	9.5mm, For .5-1.5ml Tubes	4000 Labels per Roll	Blue 
LTC-95G	9.5mm, For .5-1.5ml Tubes	4000 Labels per Roll	Green 
LTC-95N	9.5mm, For .5-1.5ml Tubes	4000 Labels per Roll	Orange 
LTC-95R	9.5mm, For .5-1.5ml Tubes	4000 Labels per Roll	Red 
LTC-95Y	9.5mm, For .5-1.5ml Tubes	4000 Labels per Roll	Yellow 
LTC-13W	13mm, For 1.5-2ml Tubes	4000 Labels per Roll	White 
LTC-13B	13mm, For 1.5-2ml Tubes	4000 Labels per Roll	Blue 
LTC-13G	13mm, For 1.5-2ml Tubes	4000 Labels per Roll	Green 
LTC-13N	13mm, For 1.5-2ml Tubes	4000 Labels per Roll	Orange 
LTC-13R	13mm, For 1.5-2ml Tubes	4000 Labels per Roll	Red 
LTC-13Y	13mm, For 1.5-2ml Tubes	4000 Labels per Roll	Yellow 


### Labels for Cryogenic Vials

Item #	Description	Quantity	Color
LTC-27X13W	27x13mm	1000 Labels per Roll	White 
LTC-27X13B	27x13mm	1000 Labels per Roll	Blue 
LTC-27X13G	27x13mm	1000 Labels per Roll	Green 
LTC-27X13N	27x13mm	1000 Labels per Roll	Orange 
LTC-27X13R	27x13mm	1000 Labels per Roll	Red 
LTC-27X13Y	27x13mm	1000 Labels per Roll	Yellow 
LTC-38X13W	38x13mm	750 Labels per Roll	White 
LTC-38X13B	38x13mm	750 Labels per Roll	Blue 
LTC-38X13G	38x13mm	750 Labels per Roll	Green 
LTC-38X13N	38x13mm	750 Labels per Roll	Orange 
LTC-38X13R	38x13mm	750 Labels per Roll	Red 
LTC-38X13Y	38x13mm	750 Labels per Roll	Yellow 

### Labels for Microplates

Item #	Description	Quantity	Color
LTC-51X6W	51x6mm	1000 Labels per Roll	White 

### Labels for Microscope Slides

Item #	Description	Quantity	Color
LTC-22X22W	22x22mm	1000 Labels per Roll	White 














# Cryogenic and Cold Storage Labels

## Direct Thermal Label Rolls











LTC-38X19W

### Labels for Large Tubes, Racks & Boxes

Item #	Description	Quantity	Color
LTC-38X19W	38x19mm	750 Labels per Roll	White 
LTC-38X19B	38x19mm	750 Labels per Roll	Blue 
LTC-38X19G	38x19mm	750 Labels per Roll	Green 
LTC-38X19N	38x19mm	750 Labels per Roll	Orange 
LTC-38X19R	38x19mm	750 Labels per Roll	Red 
LTC-38X19Y	38x19mm	750 Labels per Roll	Yellow 
LTC-43X19W	43x19mm	500 Labels per Roll	White 
LTC-51X19W	51x19mm	500 Labels per Roll	White 
LTC-51X19B	51x19mm	500 Labels per Roll	Blue 
LTC-51X19G	51x19mm	500 Labels per Roll	Green 
LTC-51X19N	51x19mm	500 Labels per Roll	Orange 
LTC-51X19R	51x19mm	500 Labels per Roll	Red 
LTC-51X19Y	51x19mm	500 Labels per Roll	Yellow 

### Label/Dot Combo for Tubes

Item #	Description	Quantity	Color
LTC-24X13-95W	24x13mm & 9.5mm Dot, For .5-1.5ml Tubes	500 Labels per Roll	White  
LTC-24X13-95B	24x13mm & 9.5mm Dot, For .5-1.5ml Tubes	500 Labels per Roll	Blue  
LTC-24X13-95G	24x13mm & 9.5mm Dot, For .5-1.5ml Tubes	500 Labels per Roll	Green  
LTC-24X13-95N	24x13mm & 9.5mm Dot, For .5-1.5ml Tubes	500 Labels per Roll	Orange  
LTC-24X13-95R	24x13mm & 9.5mm Dot, For .5-1.5ml Tubes	500 Labels per Roll	Red  
LTC-24X13-95Y	24x13mm & 9.5mm Dot, For .5-1.5ml Tubes	500 Labels per Roll	Yellow  
LTC-38X13-95W	38x13mm & 9.5mm Dot, For 1.5-2ml Tubes	500 Labels per Roll	White  
LTC-38X13-95B	38x13mm & 9.5mm Dot, For 1.5-2ml Tubes	500 Labels per Roll	Blue  
LTC-38X13-95G	38x13mm & 9.5mm Dot, For 1.5-2ml Tubes	500 Labels per Roll	Green  
LTC-38X13-95N	38x13mm & 9.5mm Dot, For 1.5-2ml Tubes	500 Labels per Roll	Orange  
LTC-38X13-95R	38x13mm & 9.5mm Dot, For 1.5-2ml Tubes	500 Labels per Roll	Red  
LTC-38X13-95Y	38x13mm & 9.5mm Dot, For 1.5-2ml Tubes	500 Labels per Roll	Yellow  
LTC-33X13-95W	33x13mm & 9.5mm Dot, For 1.5-2ml Tubes	500 Labels per Roll	White  
LTC-33X13-95B	33x13mm & 9.5mm Dot, For 1.5-2ml Tubes	500 Labels per Roll	Blue  
LTC-33X13-95G	33x13mm & 9.5mm Dot, For 1.5-2ml Tubes	500 Labels per Roll	Green  
LTC-33X13-95N	33x13mm & 9.5mm Dot, For 1.5-2ml Tubes	500 Labels per Roll	Orange  
LTC-33X13-95R	33x13mm & 9.5mm Dot, For 1.5-2ml Tubes	500 Labels per Roll	Red  
LTC-33X13-95Y	33x13mm & 9.5mm Dot, For 1.5-2ml Tubes	500 Labels per Roll	Yellow  



LTC-33X13-95W