



Task Light

LE 50[®] LED

Burton's economical multi-purpose LED light.

Burton's LE 50 LED is a cost-effective, flex-arm, medical task, and reading light. Its illumination far exceeds the competition.



burton[®]

Applications

Medical Task

- Ambulatory Surgery Centers
- Urgent Care
- Emergency / Triage Centers
- Long Term Care
- Veterinary
- School Nurse


Patient Reading

- Patient Rooms
- Senior Living / Long Term Care

Features

- High-quality illumination
- Durable construction
- Energy efficient, low maintenance LED
- Easy positioning light head
- Convenient on/off switch
- Standard or wide beam models available
- Mobile floor stand, or space-saving wall or table mount models available

LE 50® LED product dimensions



Summary of technical data

LE 50® LED

| | |
|------------------------------|---|
| Illuminance: | Standard Beam: 40,000 Lux at 18"; 75,000 Lux at 12" Wide Beam: 9,000 Lux at 18"; 19,000 Lux at 12" |
| Color Temperature: | 4100 Kelvin |
| CRI (Color Rendering Index): | 95 |
| Power Supply: | Universal Voltage 100/230V External Adapter |
| Power: | 7 Watts |
| Light Source: | Light Emitting Diode (LED) |
| Rated Life of LED: | 50,000-hour average rated LED life |
| Mounting Options: | Floor Stand, Wall, Table |
| Arm: | Flexible Gooseneck |
| Certifications/Approvals: | UL / cUL |
| Product Warranty: | 1-year limited warranty |
| Origin: | Assembled in USA |

Ordering Information

LE 50® LED - STANDARD BEAM

| Universal Input 100V/230V | Description |
|------------------------------|--|
| LE50FL | LE 50 LED Standard Beam Mobile Floor Stand |
| LE50W | LE 50 LED Standard Beam Universal Wall Mount |
| LE50T | LE 50 LED Standard Beam Table Mount |


LE 50® LED - WIDE BEAM

| | |
|---------|--|
| LEW50FL | LE 50 LED Wide Beam Mobile Floor Stand |
| LEW50W | LE 50 LED Wide Beam Universal Wall Mount |
| LEW50T | LE 50 LED Wide Beam Table Mount |

ACCESSORIES

| Model | Description |
|---------|--------------|
| 5000167 | Wall Bracket |

5000167



Mounting holes on 4 corners

© 2021 Burton Medical. All rights reserved. Burton Medical reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

burton®