

**RADIATION PROTECTIVE SHIELD
FOR STRYKER FLUORO STRETCHER**

PRODUCT INFORMATION

PRODUCT INFORMATION

RADIATION PROTECTIVE SHIELD FOR STRYKER FLUORO STRETCHER



Enhanced Radgenic Radiation Protection: X-Ray Shield for Stryker Fluoro Stretchers

Optimize patient and staff safety during fluoroscopic procedures with our custom-designed X-ray shield, specifically engineered for seamless integration with Stryker Fluoro Stretchers.

This essential accessory provides superior radiation attenuation, effectively reducing scatter radiation exposure for clinicians and nearby personnel in the OR, ER, and other imaging environments. Crafted from high-quality, .5mm lead-equivalent materials, our shield offers robust protection without impeding workflow or patient access.

Key Features & Benefits:

- Precision Fit for Stryker Fluoro Stretchers: Engineered to perfectly complement the design of Stryker Fluoro stretchers, ensuring secure attachment and optimal coverage without interfering with stretcher articulation or patient positioning.
- Superior Radiation Attenuation: Constructed with advanced lead-equivalent materials, this shield significantly reduces scatter radiation, protecting healthcare professionals from harmful exposure during lengthy or complex fluoroscopic cases.
- Optimized Workflow & Access: Thoughtfully designed to provide maximum radiation protection while maintaining critical patient access and allowing for efficient clinical procedures. Its intuitive design ensures quick and easy deployment and retraction.
- Durable & Easy to Clean: Built with Phillips exclusive Heat Sealed Non Porous Radgenic Fabric – with no sewn seam design – to effectively clean and disinfect with Cavi Wipes, promoting a safe clean environment.
- Increased Staff Confidence: By minimizing radiation exposure, the shield enhances the safety and confidence of medical staff, allowing them to focus on providing optimal patient care.
- Patient-Centric Design: While prioritizing staff safety, the shield's design also considers patient comfort and accessibility, ensuring a smooth and unencumbered experience.
- Configuration: Sold with both Side and Head/Foot shields or order separately

Invest in the safety of your team and enhance the efficiency of your fluoroscopic procedures with our Radgenic X-ray shield, the ideal companion for your Stryker Fluoro stretcher.

SPECIFICATION

MATERIAL: Radgenic®

COLOR: Blue

LEAD EQUIVALENCY: 0.50mm Pb Eq

CONFIGURATION: Sold with both Side and Head/Foot shields or order separately