

**ECONOMY L BLOCK SHIELD
CLEAR 2.0MM PB RADIATION SHIELD**

PRODUCT INFORMATION



**PHILLIPS SAFETY
PRODUCTS INC.**

PRODUCT INFORMATION

ECONOMY L BLOCK SHIELD - CLEAR 2.0MM PB RADIATION SHIELD



Economy L Block Shield – Clear-Pb Radiation Shield

Protect yourself with confidence using the Economy L Block Shield. This tabletop “L”-shaped barrier is constructed from high-quality Clear-Pb radiation shielding glass, offering 2.0mm Pb equivalency for reliable protection against low-energy gamma and beta radiation. Designed for lab and medical environments, it delivers excellent visibility while keeping your face and upper body fully shielded.

- 2.0mm Pb equivalency for dependable radiation protection
- Clear-Pb construction offers full visibility of the work area
- Portable and durable, ideal for mobile or compact workstations
- Ergonomic L-shape allows full accessibility to the workspace
- Rugged, reliable, and affordable—this shield is the perfect balance of safety and value for routine radiation protection.

SPECIFICATION

SHAPE: L-Block

FRAME SIZE: 18"W x 24"H

WINDOW SIZE: 16.5"W x 22.625"H

BASE SIZE: 14"W x 4"H

LEAD EQUIVALENCY: 2.0mm Pb Eq

FEATURES: Portable, Durable, Ergonomic L-Shape, Rugged, Reliable, Affordable

PRODUCT INFORMATION