



Technical sheet

Sterisystem® Perfo-Safe® Wire Shelves

Article number: **30003 serie**

Our Sterisystem® highly polished Perfo-Safe® wire shelves are the perfect solution if you need to sterilise, transport and store larger sterilisation packs. Manufactured in accordance with the highest quality standards and finished to perfection, these wire trays meet the most stringent CSSD and OR requirements. Our wire shelves can be used in all our transport carts and storage systems.

Sizing					
	30003 00003 DIN 1/1 (Deborah)	30003 00005 DIN 2/1 (Dorothe)	30003 00004 ISO 1/1 (Doris)	30003 00007 ISO 1/1 raised edge (Delphine)	19031 00031*
Width	23 ⁵ / ₈ inch	23 ⁵ / ₈ inch	23 ⁵ / ₈ inch	23 ⁵ / ₈ inch	950 mm
Depth	11 ¹³ / ₁₆ inch	23 ⁵ / ₈ inch	15 ³ / ₄ inch	15 ³ / ₄ inch	600 mm
Height	-	-	-	1 ⁴⁹ / ₆₄ inch	18 mm

* Only for use in the top compartment of the Sterisystem™ closed stainless steel case cart (19031 00030)

Capacity	
General capacity	55 lbs (for ergonomic reasons we advise you not to load more than 17 lbs on each shelf)

Cleaning	
Resistant to common cleaning and disinfection products used in health care. It is recommended to test any cleaning solution on a small, inconspicuous area to ensure no damage occurs as a result. For more information, please refer to the document "Belintra Cleaning Recommendations".	

Materials	
General	AISI 304L stainless steel - electrolytically polished


Temperature range	
General	Can be used inside the autoclave

Packaging information					
	30003 00003 DIN 1/1 (Deborah)	30003 00005 DIN 2/1 (Dorothe)	30003 00004 ISO 1/1 (Doris)	30003 00007 ISO 1/1 raised edge (Delphine)	19031 00031
Net/unit	2.24 lbs	3.87 lbs	2.52 lbs	2.75 lbs	12.02 lbs
Gross packaged	46.96 lbs	41.23 lbs	26.98 lbs	29.76 lbs	75.84 lbs
Units/box	20	10	10	10	6
Boxes/pallet	-	40	-	40	10


Accessories	
Silicon mat - Standard	30010 00004
Silicon mat - Heavy Duty	30010 00006

Technical drawings - see back


» 30003 00003 - DIN 1/1 (Deborah)




» 30003 00005 - DIN 2/1 (Dorothe)



» 30003 00004 - ISO 1/1 (Doris)



» 30003 00007 - ISO 1/1 raised edge (Delphine)



» 19031 00031

Only for use in the top compartment of the Sterisystem™ closed stainless steel case cart (19031 00030)

