

M7

Bladder Scanner



10,000

Cases in Data Storage

2,000ml

Max Scanning Range

128

Sensors

≥3

Seconds Scan Time

Features

- $\pm 5\%$ Accuracy
- 10.3" Color TFT-LCD Screen
- Durable, impact-resistant probes
- Patient ID support: Up to 18 digits
- Battery-powered with DC charging option
- Innovative negative pressure probe technology
- Quick guide and step by step instructions included on screen
- Search for a patient by name to quickly review all historical readings for that patient

Scan Modes: 3D Scan

Print Mode: Wireless Pinter

Working Mode(s): Easy & Expert Mode

Accessories: Medical rolling cart & travel case included. Bluetooth thermal printer optional.

Precision Scanning Technology

- Single motor mechanical & electronic hybrid scanning for efficiency
- Automatic Patient Type Selection
- **Adjustable Scanning Depths:**
160mm, 180mm, & 200mm
- **Wide Scanning Range:**
96° Swing Angle & 180° Rotation Angle



Primary Care



Longterm Care