



REVISION DATE: 3/04

I. PRODUCT IDENTIFICATION: OB TOWELETTE - D74800

PRODUCT/TRADE NAME:

HAZARD RATING(NFPA)

HEALTH: 0

FLAMMABILITY: 1

REACTIVITY: 0

SPECIFIC: NONE

EMERGENCY OR INFORMATION TELEPHONE NO: 845-365-1700 (M-F DAYTIME) AT OTHER TIMES, CONTACT THE LOCAL POISON CONTROL CENTER

CHEMICAL NAME: PROPRIETARY MIXTURE

II HAZARDOUS INGREDIENTS PER 29 CFR 1910.1200

Hazardous Ingredients	%	ACGIH TLV	CAS NUMBER
SD ALCOHOL 40	5	1000 PPM	64-17-5

III-PHYSICAL/CHEMICAL CHARACTERISTICS

III-PHYSICAL/CHEMICAL CHARACTERISTICS

COLOR/ODOR/APPEARANCE: TOWELETTE SATURATED WITH

HAZY LIQUID WITH THYMOL ODOR..

BOILING POINT: N/A

FLASH POINT: 169°F

VAPOR DENSITY N/A

EVAPORATION RATE: N/A

SOLUBILITY IN WATER: COMPLETE

SPECIFIC GRAVITY (H₂O = 1): 0.994

IV FIRE & EXPLOSION HAZARD DATA

FLASH POINT(Method Used): 169°F TAG CLOSED CUP : LEL: N/A UEL: N/A

EXTINGUISHING MEDIA: DRY CHEMICAL OR ALCOHOL TYPE FOAM, CARBON DIOXIDE

SPECIAL FIRE FIGHTING PROCEDURES: NONE

USUAL FIRE AND EXPLOSION HAZARDS: NONE

V-REACTIVITY DATA

STABILITY: STABLE

CONDITIONS TO AVOID NONE

INCOMPATIBILITY: NONE

HAZARDOUS DECOMPOSITION OR BYPRODUCT: NONE

POLYMERIZATION: WILL NOT OCCUR

CONDITIONS TO AVOID: NONE

OB TOWELETTE

VI-HEALTH HAZARD DATA

EFFECTS OF OVER EXPOSURE:

SKIN IF RASH APPEARS, DISCONTINUE USE

EYES : WILL CAUSE EYE STING IF SPLASHED

INHALATION: NONE

INGESTION NOT A NORMAL ROUTE OF ENTRY

EMERGENCY AND FIRST AID PROCEDURES:

SKIN CONTACT: DISCONTINUE USE IF RASH OR IRRITATION OCCURS

EYE CONTACT: FLUSH WITH COLD WATER FOR 15 MINUTES

INHALATION: REMOVE TO FRESH AIR

INGESTION: IF SWALLOWED, CALL POISON CONTROL CENTER

TARGET ORGANS; N/A

VII-SPILL AND DISPOSAL PROCEDURE

SPILL CONTROL: ELIMINATE ALL SOURCES OF IGNITION

WASTE DISPOSAL METHOD: FLUSH SPILLS WITH WATER. HANDLE AS FLAMMABLE LIQUID. FOLLOW LOCAL, STATE AND FEDERAL REGULATIONS.

HANDLING AND STORAGE: STORE AWAY FROM HEAT AND SOURCES OF IGNITION

VIII-CONTROL MEASURES/PROTECTION

RESPIRATION: USE IF PERMISSABLE EXPOSURE LEVEL IS EXCEEDED WHEN HANDLING BULK LIQUID

VENTILATION: RECOMMENDED

PROTECTIVE GLOVES: NONE REQUIRED

EYE PROTECTION: ONLY IF SPLASHING IS EXPECTED

HYGIENIC PRACTICES:GOOD HOUSEKEEPING PRACTICES SHOULD BE FOLLOWED.

OTHER: DON'T ALLOW LARGE QUANTITIES OF WASTE TO ACCUMULATE

IX TRANSPORT/SHIPPING

DOT SHIPPING NAME: CONSUMER COMMODITY

TECHNICAL SHIPPING NAME: N/A

DOT SHIPPING CLASSIFICATION: DOT NOT REGULATED

DOT ID NO.: N/A

DOT LABEL REQUIREMENTS: N/A

UN/NA NUMBER REGULATIONS: N/A

REPORTABLE QUANTITY: N/A

X DISCLAIMER

THE INFORMATION FURNISHED HEREIN IS BELIEVED TO BE ACCURATE AND REPRESENTS THE BEST DATA CURRENTLY AVAILABLE TO US. NO WARRANTY, EXPRESSED OR IMPLIED IS MADE AND NICE-PAK PRODUCTS, INC. ASSUMES NO LEGAL RESPONSIBILITY OR LIABILITY RESULTING FROM ITS USE.