

# EZVISION SERIES

## Video Laryngoscope

---

A fully portable, high-performance solution for emergency care, designed with ergonomic excellence in mind! Seamlessly connect to larger monitors and computers, allowing for real-time photo and video capture. With its compact size and remarkable versatility, the EZVision series becomes an indispensable tool for emergency care providers.



**.55 lbs**

Light  
Weight

**5 HR**

Continuous  
Usage

**3"**

HD LCD  
Display

# Features

---

- Clear visualization
- 3" Wide angle LCD video
- Portable, lightweight, & waterproof (IPX7)
- Anti-fog blades with locking feature
- Disposable blades preventing cross contamination
- Real-time capture allows for instant capture of photos and videos
- Versatile Connectivity: Can be effortlessly connected to larger monitors and computers

## **EZVision Blades:**

Pediatric MAC 2, Medium Adult MAC 3, Large Adult MAC 4

## **EZVision+ Blades:**

Premature Mil 00, Neonate Mil 0, Pediatric Mil 1

# Ergonomic Design

---

Comfortable grip and intuitive controls, allowing for **precise** and **efficient airway management**.

# Market Type

---



Acute



ASC



EMS



Primary Care



**EZVision Plus**  
(See Page 4)

## EZVision Blades

---



**Pediatric MAC 2**  
*Part Number: YLO-III B2*




**Medium Adult MAC 3**  
*Part Number: YLO-III B3*



**Large Adult MAC 4**  
*Part Number: YLO-III B4*

## Sample Images

---



## Specifications

---

### Display

Size: 3"  
Resolution: 1280 x 1024px  
Power: 3.7V

### Laryngoscope Blade (Including Camera)

Resolution: 2 mega-pixels  
Power: 3.7V  
Light Source: LED 800LUX

### Battery

Type: Rechargeable lithium battery  
Charging times: >500 times  
Charging time: <5 hours  
Operating time: Lasting 5 Hours

## Blades & Stands

---




**MIL 00 (Premature)**  
Part Number: YLO-IB00




**MIL 0 (Neonate)**  
Part Number: YLO-IB0



**MIL 1 (Child)**  
Part Number: YLO-IB1



**MIL 00/0 Stand (Premature/Neonate)**  
Part Number: YLO-IS0  
*\*INCLUDED IN PURCHASE\**



**MIL 1 Stand (Pediatrics)**  
Part Number: YLO-IS1  
*\*INCLUDED IN PURCHASE\**

## EZVision+

---

The EZVision+ is a video laryngoscope for premature and neonate demographics that has interchangeable stands for Miller 00, 0, & 1 sizings. Features disposable anti-fog blades for clear visibility and 3" LCD Screen.



**Detachable LCD Screen**  
Part Number: EZVision Plus

## Specifications

---

### Display

Size: 3"  
Resolution: 1600 x 1200px  
Power: 3.7V

### Laryngoscope Blade

Luminance: LED  $\geq 150$  LUX  
Pixels: 2 Mega-Pixels  
Power: 3.7V

### Battery

Power: 3.7V  
Type: Rechargeable lithium battery  
Charging times: >500 times  
Charging time: <4 hours  
Operating time: Lasting 5 Hours