

MATERIAL SAFETY DATA SHEET

PRODUCT IDENTIFICATION

PRODUCT/TRADE NAME: SANI CLOTH HB Q10584, **Q85484**, Q08472, H21082, U22495HAZARD RATING(NFPA)

HEALTH: 0
FLAMMABILITY: 0
REACTIVITY: 0
SPECIFIC: NONE

EMERGENCY OR INFORMATION TELEPHONE NO: 845-365-1700 (M-F DAYTIME) AT OTHER TIMES, CONTACT THE LOCAL POISON CONTROL CENTER

CHEMICAL NAME: PROPRIETARY MIXTURE

II HAZARDOUS INGREDIENTS PER 29 CFR 1910.1200

<u>Hazardous Ingredients</u>	<u>%</u>	<u>ACGIH TLV</u>	<u>CAS NUMBER</u>
QUATERNARY AMMONIUM CMPD	1.0	N/A	68391-01-5/68956-79-6

III PHYSICAL/CHEMICAL CHARACTERISTICS

COLOR/ODOR/APPEARANCE: APPLICATOR SATURATED WITH CLEAR COLORLESS LIQUID WITH MILD QUAT ODOR
BOILING POINT: N/A
FLASH POINT: NOT APPLICABLE
VAPOR DENSITY: N/A
EVAPORATION RATE: N/A
SOLUBILITY IN WATER: COMPLETE
SPECIFIC GRAVITY (H₂O = 1): 1.0

IV FIRE & EXPLOSION HAZARD DATA

FLASH POINT(Method Used): NOT APPLICABLE : LEL: N/A UEL: N/A
EXTINGUISHING MEDIA: NONE
SPECIAL FIRE FIGHTING PROCEDURES: NONE
UNUSUAL FIRE AND EXPLOSION HAZARDS: NONE

V-REACTIVITY DATA

STABILITY: STABLE
CONDITIONS TO AVOID : NONE
INCOMPATIBILITY: ANIONIC SURFACTANT, SOAPS
HAZARDOUS DECOMPOSITION OR BYPRODUCTS: NONE KNOWN
POLYMERIZATION: WILL NOT OCCUR
CONDITIONS TO AVOID: NONE KNOWN

VI-HEALTH HAZARD DATA

EFFECTS OF OVER EXPOSURE:SKIN : MAY CAUSE SKIN IRRITATION
EYES : MAY CAUSE IRRITATION
INHALATION: NOT NORMAL ROUTE OF ENTRY
INGESTION: NOT NORMAL ROUTE OF ENTRY

EMERGENCY AND FIRST AID PROCEDURES:

SKIN CONTACT: IF RASH OR IRRITATION DEVELOPS, DISCONTINUE USE
EYE CONTACT: FLUSH WITH COOL WATER FOR 15 MINUTES
INHALATION: REMOVE TO FRESH AIR
INGESTION: CONSULT PHYSICIAN
TARGET ORGANS; NONE KNOWNKNOWN
CARCINOGEN: NONE KNOWN

VII-SPILL AND DISPOSAL PROCEDURE

SPILL CONTROL: PICK UP WITH ABSORBENT
WASTE DISPOSAL METHOD: OBEY LOCAL, STATE AND FEDERAL REGULATIONS
HANDLING AND STORAGE: STORE AT ROOM TEMPERATURE

VIII-CONTROL MEASURES/PROTECTION

RESPIRATION: USE NIOSH APPROVED RESPIRATOR IF EXPOSURE LEVEL IS EXCEEDED
VENTILATION: RECOMMENDED
PROTECTIVE GLOVES: RECOMMENDED
EYE PROTECTION: IF SPLASHING MAY OCCUR
HYGIENIC PRACTICES: FOLLOW GOOD HOUSEKEEPING PRACTICES
OTHER: NONE

IX TRANSPORT/SHIPPING

DOT SHIPPING NAME:
TECHNICAL SHIPPING NAME:
DOT SHIPPING CLASSIFICATION: DOT NOT REGULATED
DOT ID NO.:
DOT ;LABEL REQUIREMENTS:NA
UN/NA NUMBER REGULATIONS:N/A
REPORTABLE QUANTITY: N/A

X DISCLAIMER

THE INFORMATION FURNISHED HEREIN IS BELIEVED TO BE ACCURATE AND REPRESENTS THE BEST DATA CURRENTLY AVAILABLE TO US. NO WARRANTY, EXPRESSED OR IMPLIED IS MADE AND NICE-PAK PRODUCTS, INC. ASSUMES NO LEGAL RESPONSIBILITY OR LIABILITY RESULTING FROM ITS USE.