



*fixed angle centrifugation

DASH Centrifuges - High Performance STAT Processing

Model	Rotor Type	Interface	Max Speed RPM / RCF	Up to 13 x 75 mm or 5 mL	Up to 16 x 100 mm or 10 mL	Up to 16 x 125 mm or 15 mL	50 mL	Part Number
DASH Apex 4	Fixed Angle	2-Button, 3-Setting	8500 / 4440	4	-	-	-	00-053-009-000
DASH Coag	Fixed Angle	2-Button, 3-Setting	6600 / 4400	-	12	-	-	00-676-009-001
DASH Flex 6	Horizontal	Digital, 10-Setting	5400 / 4150	-	6	-	-	00-176-009-000
DASH Flex 12	Horizontal	Digital, 10-Setting	5200 / 4000	-	12	-	-	00-183-009-000
DASH Flex 24	Horizontal	Digital, 10-Setting	4500 / 3450	-	24	12	6	00-184-009-000

(adapters available for a variety of tube sizes)

SERO 12 - Blood Banking Centrifuge A drop-in replacement for the BD SeroFuge

Model	Rotor Type	Interface	Max Speed RPM / RCF	12 Place Rotor Up to 75mm	6 Place Rotor Up to 100 mm	Part Number
SERO 12	Fixed Angle	Digital, 10-Setting	3600 / 1040	6 capped 12 uncapped	6	00-483-009-001

HORIZON Centrifuges Quality Routine Chemistry

Model	Rotor Type	Interface	Max Speed RPM / RCF	Up to 16 x 100 mm or 10 mL	Up to 16 x 125 mm or 15 mL	50 mL	Part Number
HORIZON 6 Flex	Horizontal	Digital, 10-Setting	3800 / 2000	6	-	-	00-376-009-000
HORIZON 6 Flex FA	Fixed Angle	Digital, 10-Setting	3900 / 1850	6	6	-	00-376-009-001
HORIZON 12 Flex	Horizontal	Digital, 10-Setting	3700 / 2000	12	-	-	00-383-009-000
HORIZON 24 Flex	Horizontal	Digital, 10-Setting	3400 / 2000	24	12	6	00-384-009-000

(adapters available for a variety of tube sizes)

600 Series Centrifuges Compact Routine Chemistry

Model	Rotor Type	Interface	Max Speed RPM / RCF	Up to 16 x 100 mm or 10 mL	Up to 16 x 125 mm or 15 mL	50 mL	Part Number
614B	Fixed Angle	Dial, 1-Setting	3350 / 1300	6	6	-	00-004-009-001
642B	Horizontal	Dial, 1-Setting	3380 / 1600	6	-	-	00-016-009-001
642E	Horizontal	2-Button, 1-Setting	3380 / 1600	6	-	-	00-018-009-001
642M	Horizontal	2-Button, 1-Setting	3400 / 1650	6	-	-	00-077-009-000

(adapters available for a variety of tube sizes)