



Solutions to Elevate Safety & Productivity in Healthcare.

Lead Shield Accessory Features and Benefits



Features:

- Tip to roll for easy transport
- Wheels disengage when upright for added stability
- Shield pivots 180 degrees for versatility
- Height adjustable from 49" to 16"
- 34 height adjustment options in approx. 1 inch increments
- Shield arm detaches/stows when not in use
- 1" leg profile
- Limited Lifetime warranty on materials and workmanship
- Custom shield size and color available upon request

Benefits

- Ideal for shielding abdomen and pelvis for standing x-ray exams
- Locking height adjustment reliably holds desired height
- Will not roll away during use
- 1" leg profile reduces tripping hazard
- Multi-use accessory when purchased with Leg Lift positioning device

Specifications:

Shield Dimensions: 24" wide x 23" tall (custom heights available upon request)

Weight: 9.5 lbs.

Footprint: 20" x 20" (shield extends beyond footprint)

Shipping:

Ships in one box

Ships Assembled