



Solutions to Elevate Safety & Productivity in Healthcare.

Shielded Extremity Table Accessory Features and Benefits



Features:

- Shielded table prevents x-rays from passing through
- Leg Lift base (sold separately) tips back to roll to your patient
- Wheels disengage when upright for added stability
- Height adjustable from 49" to 16"
- Extremity Table detaches for storage when not in use
- 1" leg profile minimizes trip hazard
- Limited Lifetime warranty on materials and workmanship
- 17" x 17" working surface accommodates most image receptors

Benefits

- Ideal for pediatric exams
- Can be used for seated or standing exams
- Small footprint allows for use even in small spaces
- Will not roll away during use
- Low profile feet reduce tripping hazard
- Multi-use accessory when purchased with Leg Lift positioning device

Specifications:

- Accessory depth 26.75"
- Accessory width 17.75"
- Accessory height 4.75"

Shipping:

Ships in one box

Ships assembled