



Solutions to Elevate Safety & Productivity in Healthcare.

Cassette Holder Accessory Features and Benefits



Features:

- Leg Lift base (sold separately) tips back to roll to your patient
- Wheels disengage when upright for added stability
- Cassette pivots 180 degrees for versatility
- Height adjustable from 49" to 16"
- Cassette holder detaches for storage when not in use
- 1" leg profile minimizes trip hazard
- Limited lifetime warranty on materials and workmanship

Benefits

- Small footprint allows for use even in small spaces
- Will not roll away during use
- Low profile feet reduce tripping hazard
- Multi-use accessory when purchased with Leg Lift positioning device

Specifications:

- Accessory depth 20"
- Accessory width 6"
- Accessory height 21.5"

Shipping:

Ships in one box
Ships assembled