



*Solutions to Elevate Safety & Productivity in Healthcare.*

## Leg Lift Positioning Device 2.0

### Features and Benefits

---

#### Features:

- Leg Lift arm supports up to 200 lbs.
- Tips back to roll to bedside
- Wheels disengage when upright for added stability
- Support cradle shaped for comfort
- Release clutch knob to adjust cradle for patient comfort
- 34 height adjustment options in approx. 1 inch increments
- Beveled edges to prevent skin injury
- Leg support cradle can be detached for storage
- Leg support cradle adjusts from 16" to 50"
- 1" leg profile compatible with most hospital beds and gurneys
- 1" leg profile minimizes trip hazard
- Limited Lifetime warranty on materials and workmanship

#### Benefits

- Stable platform supports 200 lbs. even for unpredictable patients
- Space beneath leg support cradle remains clear of obstruction for x-ray beam, instrumentation, or catchment system
- Reliable support for patients who cannot participate in positioning
- Moulded leg support cradle comfortably supports all leg sizes
- Indicated for cross-table lateral hip or femur projections as well as lateral and tangential knee projections, and foot and ankle exams for wheelchair bound patients. Also indicated for the operative setting, emergency departments, and other settings as approved by Flow State Innovation.

#### Specifications:

Height: 51 inches

Weight: 27 lbs

Footprint: 20" x 20"

Leg Support Cradle Height Adjustment: 16" to 50"

#### Shipping:

Ships in one 31 lbs box





*Solutions to Elevate Safety & Productivity in Healthcare.*



## Lead Shield Accessory Features and Benefits

---

### Features:

- Tip to roll for easy transport
- Wheels disengage when upright for added stability
- Shield pivots 180 degrees for versatility
- Height adjustable from 49" to 16"
- 34 height adjustment options in approx. 1 inch increments
- Shield arm detaches/stows when not in use
- 1" leg profile
- Limited Lifetime warranty on materials and workmanship
- Custom shield size and color available upon request

### Benefits

- Ideal for shielding abdomen and pelvis for standing x-ray exams
- Locking height adjustment reliably holds desired height
- Will not roll away during use
- 1" leg profile reduces tripping hazard
- Multi-use accessory when purchased with Leg Lift positioning device

### Specifications:

Shield Dimensions: 24" wide x 23" tall (custom heights available upon request)

Weight: 9.5 lbs.

Footprint: 20" x 20" (shield extends beyond footprint)

### Shipping:

Ships in one box

Ships Assembled



*Solutions to Elevate Safety & Productivity in Healthcare.*

## Shielded Extremity Table Accessory Features and Benefits

---



### Features:

- Shielded table prevents x-rays from passing through
- Leg Lift base (sold separately) tips back to roll to your patient
- Wheels disengage when upright for added stability
- Height adjustable from 49" to 16"
- Extremity Table detaches for storage when not in use
- 1" leg profile minimizes trip hazard
- Limited Lifetime warranty on materials and workmanship
- 17" x 17" working surface accommodates most image receptors

### Benefits

- Ideal for pediatric exams
- Can be used for seated or standing exams
- Small footprint allows for use even in small spaces
- Will not roll away during use
- Low profile feet reduce tripping hazard
- Multi-use accessory when purchased with Leg Lift positioning device

### Specifications:

- Accessory depth 26.75"
- Accessory width 17.75"
- Accessory height 4.75"

### Shipping:

Ships in one box

Ships assembled



*Solutions to Elevate Safety & Productivity in Healthcare.*

## Cassette Holder Accessory Features and Benefits

---



### Features:

- Leg Lift base (sold separately) tips back to roll to your patient
- Wheels disengage when upright for added stability
- Cassette pivots 180 degrees for versatility
- Height adjustable from 49" to 16"
- Cassette holder detaches for storage when not in use
- 1" leg profile minimizes trip hazard
- Limited lifetime warranty on materials and workmanship

### Benefits

- Small footprint allows for use even in small spaces
- Will not roll away during use
- Low profile feet reduce tripping hazard
- Multi-use accessory when purchased with Leg Lift positioning device

### Specifications:

- Accessory depth 20"
- Accessory width 6"
- Accessory height 21.5"

### Shipping:

Ships in one box  
Ships assembled