



A centrifuge for every animal health need.


TrueBond Centrifuges
Multipurpose Veterinary Processing

| Model | Rotor Type | Interface | Preset Cycle Time / Speed / RCF | Max Speed RPM / RCF | Up to 16 x 100 mm or 10 mL | Up to 16 x 125 mm or 15 mL | Microhematocrit | Dimensions (W x D x H) |
|----------------------------------|-------------|---------------------|---|---------------------|----------------------------|----------------------------|-----------------|---|
| TrueBond Set 3 00-276-009-007 | Fixed Angle | 2-Button, 3-Setting | 10 min / 3900 RPM / 1650 xg 5 min / 2000 RPM / 425 xg 6 min / 3800 RPM / 1550 xg | 3900 / 1850 | 6 | 6 | 8 | 12 x 14 x 9 in 30.5 x 35.6 x 22.9 cm |
| TrueBond 00-376-009-005 | Fixed Angle | Digital, 10-Setting | 10 min / 3900 RPM / 1650 xg 5 min / 2000 RPM / 500 xg 5 min / 1600 RPM / 300 xg 6 min / 3800 RPM / 1550 xg | 3900 / 1850 | 6 | 6 | 8 | 12 x 14 x 9 in 30.5 x 35.6 x 22.9 cm |

(adapters available for a variety of tube sizes)

TrueBond Duet
Compact Dual Rotor Centrifuge


| Model | Rotor Type | Interface | Preset Cycle Time / Speed / RCF | Max Speed RPM / RCF | 1.5 - 2 mL Fixed Angle Rotor | Microhematocrit HCT Rotor | Dimensions (W x D x H) |
|---------------------------------|------------|---------------------|--|--|------------------------------|---------------------------|---------------------------------------|
| TrueBond Duet 00-056-009-000 | Dual Rotor | 2-Button, 5-Setting | 30 sec / 15800 RPM / 12000 xg 1.5 min / 15800 RPM / 12000 xg 2.5 min / 15800 RPM / 12000 xg 45 sec / 9800 RPM / 3900 xg 2 min / 15800 RPM / 13700 xg | 15800 / 13700 (microhematocrit) 15800 / 12000 (1.5mL - 2mL) | 2 | 12 | 7 x 7 x 7 in 17.8 x 17.8 x 17.8 cm |

TrueBond TriFLEX
Triac Replacement Centrifuge


| Model | Rotor Type | Interface | Preset Cycle Time / Speed / RCF | Max Speed RPM / RCF | Up to 16 x 100 mm or 10 mL | Up to 16 x 125 mm or 15 mL | Microhematocrit | Dimensions (W x D x H) |
|---------------------------|-------------|---------------------|--|---------------------|----------------------------|----------------------------|-----------------|--|
| TriFLEX In Development | Combination | Digital, 10-Setting | 5 min / 3500 RPM / 1500 xg 6 min / 2300 RPM / 600 xg 5 min / 1571 RPM / 300 xg 3 min / 10400 RPM / 12600 xg | 10400 / 12600 | 8 Fixed Angle | 8 Fixed Angle | 12 Horizontal | 15 x 18 x 11 in 38.1 x 45.7 x 27.9 cm |

(adapters available for a variety of tube sizes)

