

X RAY DOSIMETER BADGE HOLDER

PRODUCT INFORMATION



PRODUCT INFORMATION

X RAY DOSIMETER BADGE HOLDER | RA-DSBG-LF50



This nylon dosimeter bag is designed to securely hold and clip your dosimeter for accurate radiation monitoring. It is ideal for professionals across medical, dental, veterinary, chiropractic, and imaging practices where quarterly or monthly dosimetry badge services are required. TSA officials, law enforcement, and security personnel who may be exposed to low levels of radiation can also benefit from its use.

Pain management clinics, hospitals, and imaging centers often rely on monthly monitoring, while smaller practices may choose quarterly service. In line with Nuclear Regulatory Commission (NRC) guidelines, any employee receiving 10% or more of the annual occupational dose limit must wear a dosimeter. Our nylon dosimeter bag ensures the badge is worn properly and consistently, providing reliable exposure tracking in compliance with safety standards.

Proudly made in the USA.

SPECIFICATION

PRODUCT TYPE: Dosimeter Bag

MATERIAL: Nylon

COLOR: Yellow

DIMENSIONS: 4" W x 5" L