


RADIATION SHIELDING PLEATED TABLE SHIELD

PRODUCT INFORMATION



PHILLIPS SAFETY
PRODUCTS INC.

PRODUCT INFORMATION
RADIATION SHIELDING PLEATED TABLE SHIELD



The Radiation Shielding Pleated Table Shield is a radiation table shield that complies with CDC Efficacy guidelines for Category IB disinfection protocol for noncritical medical devices. Available in two sizes, this Radgenic shield adds an extra layer of protection with a lead core of 0.50mm Pb Eq.

Plus, it is made from a polyurethane-coated fabric, which can be easily cleaned and disinfected with Quaternary or Isopropyl disinfectant wipes. The Radiation Shielding Pleated Table Shield is non-porous, with a smooth surface and fully sealed seams. It is also latex-free, has no Velcro, and is soil and stain-resistant. This Radgenic radiation table shield is available in blue.

SPECIFICATION

DIMENSIONS: 21.5" W x 19" L | 24" W x 24" L

LEAD EQUIVALENCY: 0.50mm Pb Eq

MATERIAL: Polyurethane coated fabric

MRI SAFE: No

COLOR: Blue

FEATURES: Non-porous, Smooth surface,
Latex free, Velcro free, Stain resistant