

Classic-Value, Cabinet Table with 4 Doors



MODEL 9019

- Same cost for adjustable backrest, or flat top (Adjustable backrest is standard. Choose 088 option for flat top)
- All laminate base construction is easier to clean and sanitize
- Ultra-durable, PVC-laminated, solid, hardwood legs
- Easy assembly with pre-threaded inserts and hanger bolts in the legs
- Two, hinged doors on each side
- Backrest has 6 settings of approximately: 60°, 55°, 50°, 43°, 36° & 27°
- Seamless rounded corner top
- 2" firm foam padding (5 cm)
- Heavyweight, knit-backed upholstery
- 400 lbs. load capacity under normal use (181 kg)
- Saves on shipping—fast, easy on-site assembly required



SPECIFICATIONS


9019-24	72" L x 24" W x 31" H
9019-27	72" L x 27" W x 31" H
9019-30	72" L x 30" W x 31" H
9019-38	78" L x 30" W x 31" H

Prop 65 Warning


WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1


Options




Laminate Finishes



Maple



Dark Cherry



Gray

Standard Upholstery Colors



Slate Blue

Wedgewood

Royal Blue

Desert Tan

Burgundy

Country Mist

Warm Gray

Gunmetal

Black