

Wedge



MODEL 56

- Firm high density urethane foam core
- 18° incline
- Knit-backed, sewn, vinyl cover with hook and loop closure
- Easy-clean surface


NOTE: If no vinyl color is selected, Royal Blue will be used.

MODEL 56 SPECIFICATIONS

24" W x 18" D x 6" H at tall end



California residents please be advised of the following, pursuant to Proposition 65:

 **WARNING:** This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to www.P65Warnings.ca.gov.

2

Standard Upholstery



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray




Gunmetal



Black

Wedge



MODEL 57

- Firm high density urethane foam core
- 20° incline
- Knit-backed, sewn, vinyl cover with hook and loop closure
- Easy-clean surface

NOTE: If no vinyl color is selected, Royal Blue will be used.

MODEL 57 SPECIFICATIONS

20" W x 32" D x 12" H at tall end



California residents please be advised of the following, pursuant to Proposition 65:

! **WARNING:** This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to www.P65Warnings.ca.gov.

2

Standard Upholstery



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black

Wedge



MODEL 58

- Firm high density urethane foam core
- 65° slopes
- Knit-backed, sewn, vinyl cover with hook and loop closure
- Easy-clean surface

NOTE: If no vinyl color is selected, Royal Blue will be used.

MODEL 58 SPECIFICATIONS

16" L x 8" W x 8" H



California residents please be advised of the following, pursuant to Proposition 65:

! **WARNING:** This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to www.P65Warnings.ca.gov.

2

Standard Upholstery

