

## Cabinet Table with Center Drawers & Open Storage



### MODEL 98013

- Choose standard adjustable backrest or flat top option at no additional charge
- All laminated construction
- Backrest has 6 settings of approximately: 60°, 55°, 50°, 43°, 36° & 27°
- Center drawer section has three drawers that hold up to 75 lbs.
- Drawer section is easily reversible for access from either side
- Drawers have smooth glides for easy open and close
- Open storage with shelf on each side of center section
- Seamless rounded corner top
- 2" firm foam padding
- Heavyweight, knit-backed upholstery
- Ships assembled
- 425 lbs. load capacity under normal use (193 kg)

### SPECIFICATIONS

98013-27	72" L x 27" W x 31" H
98013-30	72" L x 30" W x 31" H
98013-28	78" L x 27" W x 31" H
98013-38	78" L x 30" W x 31" H

### Prop 65 Warning

**!** **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

### Options



### Laminate Colors



### Upholstery Colors

