

Cabinet Table with 4 Sliding Doors & Open Storage



MODEL 98044


- Choose standard adjustable backrest or flat top option at no additional charge
- All laminated construction
- Backrest has 6 settings of approximately: 60°, 55°, 50°, 43°, 36° & 27°
- Ultra-durable PVC edges
- Sliding doors can be positioned for open shelf access left, right or center
- Sliding doors are on both sides of the table
- Each section has one shelf
- Easy-access with doors on both sides of the table
- Seamless rounded corner top
- 2" firm foam padding
- Heavyweight, knit-backed upholstery
- Ships assembled
- 425 lbs. load capacity under normal use (193 kg)



SPECIFICATIONS

98044-27	72" L x 27" W x 31" H
98044-30	72" L x 30" W x 31" H
98044-28	78" L x 27" W x 31" H
98044-38	78" L x 30" W x 31" H

Prop 65 Warning



WARNING: This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to www.P65Warnings.ca.gov.

2

Options



Laminate Colors



Upholstery Colors

