

Cabinet Table with Center Open Shelf & Closed Door Storage



MODEL 98074

- Choose standard adjustable backrest or flat top option at no additional charge
- All laminated construction
- Backrest has 6 settings of approximately: 60°, 55°, 50°, 43°, 36° & 27°
- Two, closed-door, storage compartments with access from both sides
- 1 adjustable shelf behind each door
- Center open storage with shelf
- Seamless rounded corner top
- 2" firm foam padding
- Heavyweight, knit-backed upholstery
- Ships assembled
- 425 lbs. load capacity under normal use (193 kg)

SPECIFICATIONS

98074-27	72" L x 27" W x 31" H
98074-30	72" L x 30" W x 31" H
98074-28	78" L x 27" W x 31" H
98074-38	78" L x 30" W x 31" H

Prop 65 Warning

WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1

Options



20- Small Pillow

60- Large Pillow

040- Paper Cutter

030- Paper Dispenser

034- Paper Dispenser & Cutter Combo

009- Nose Hole

055- Door Locks

088- Flat Top

Laminate Colors



Maple



Gray

Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black