

# VEST AND SKIRT MOBILE RACK

PRODUCT INFORMATION

## PRODUCT INFORMATION

### VEST AND SKIRT MOBILE RACK



Efficiently store and transport your radiation skirt and vest aprons with the Phillips Safety 5 Apron Mobile Rack. Designed for medical and imaging environments, this rack securely holds up to five complete skirt and vest aprons—keeping them organized, accessible, and free from creasing or damage.

The upper angled arms are designed to support vests, while the lower horizontal rods accommodate matching skirts. Each arm swings in both directions for easy loading and removal. Built from heavy-duty materials and equipped with durable locking casters, the rack moves smoothly across any floor surface while safely supporting the full weight of radiation protection garments.

#### Specifications:

- Capacity: Holds 5 complete skirt and vest aprons
- Dimensions: 25.5 L x 25.5" W x 72" H
- Weight: 78 lbs
- Arm Length: 23"
- Swing arms move in both directions for convenient access

This mobile apron rack is the ideal solution for hospitals, radiology departments, and imaging centers seeking to protect and extend the lifespan of their valuable radiation protection garments.

## SPECIFICATION

**BRAND:** Phillips Safety

**PRODUCT TYPE:** Lead Apron Rack

**RACK TYPE:** Mobile Rack

**CAPACITY:** 5 Vest and Skirt Aprons

**WEIGHT:** 78 lbs

**ARM LENGTH:** 23"

**OVERALL SIZE:** 25.5 L x 25.5" W x 72" H

**FEATURES:** Upper angled arms for vests; Lower horizontal rods for skirts; Durable locking casters; Arm swings in both directions

## PRODUCT INFORMATION

**PRODUCT INFORMATION**  
**VEST AND SKIRT MOBILE RACK**

