

# AM-BA SERIES

## 34" BARIATRIC HI-LO TREATMENT TABLES

### STANDARD FEATURES:

- Load capacity: 500 lbs
- 34"W x 76"L x 18-37"H
- 34 oz. Heavy-duty vinyl with Permablok 3® bacterial protection
- 1.5" High density foam
- Rounded corners
- Steel base with dual foot bar control
- White powder coat base
- (4) Total locking casters
- ADA Compliant



AM-BA234

### TABLE CONFIGURATIONS:

<input type="radio"/> AM-BA134 - 1 Section	<input type="radio"/> AM-BA234 - 2 Section	<input type="radio"/> AM-BA334 - 3 Section
QTY	QTY	QTY



### TABLE FEATURES:

<p>Flat Cushion</p>	<p>Lower Extremity Angle Ranges from 0 to 75°</p>	<p>Lower Extremity Angle Ranges from 0 to 78°</p> <p>Head Section Angle Ranges from +45° to -90°</p>
AM-BA134	AM-BA234	AM-BA334

**OPTIONS:**

<input type="radio"/> <b>AM-17</b> - Face Plug
<input type="radio"/> <b>AM-91</b> - Add 1" Extra Soft Foam
<input type="radio"/> <b>19859</b> - Add-on Belt & Brackets
<input type="radio"/> <b>SPCO</b> - Special Vinyl Color (Provided by Customer)
<input type="radio"/> <b>SPBASE</b> - Special Base Color (Select Custom Base Color)
<input type="radio"/> <b>12088</b> - Side Panel 9 x 33.25 Perforated
<input type="radio"/> <b>12075</b> - Side Panel 9 x 33.25 Flat

**STANDARD UPHOLSTERY CHOICES:**

<input type="radio"/> Imperial Blue	<input type="radio"/> Taupe
<input type="radio"/> Blue Ridge	<input type="radio"/> Burgundy
<input type="radio"/> Forest Green	<input type="radio"/> Dove Gray
<input type="radio"/> Black	<input type="radio"/> Greystone
<input type="radio"/> Cappuccino	<input type="radio"/> Merlot
<input type="radio"/> Patina	

**STANDARD BASE COLOR:**

<input type="radio"/>	Office White
-----------------------	--------------

**CUSTOM BASE COLORS:**

<input type="radio"/>	Black
<input type="radio"/>	Steel Gray
<input type="radio"/>	Metallic Silver