

AM SERIES

34" OR 36" BARIATRIC HI-LO TREATMENT TABLES

STANDARD FEATURES:

- Load capacity: 500 lbs (AM-134, AM-234, AM-334)
- Load capacity: 800 lbs (AM-368 & AM-368PB)
- 34" or 36" W x 76" L x 18-37" H
- 34 oz. Heavy-duty vinyl with Permablok 3® bacterial protection
- 1.5" High density foam
- Rounded corners
- Foot control pedal (hand control on AM-368)
- White powder coat base
- ADA Compliant



AM-368

TABLE CONFIGURATIONS:

<input type="radio"/> AM-134 - 1 Section	<input type="radio"/> AM-234 - 2 Section	<input type="radio"/> AM-334 - 3 Section	<input type="radio"/> AM-368 - 3 Section Super Bariatric	<input type="radio"/> AM-368PB - 3 Section Super Bariatric with Power Tilt Back
QTY	QTY	QTY	QTY	QTY














TABLE FEATURES:

<p>Flat Cushion</p>	<p>Lower Extremity Angle Ranges from 0 to 75°</p>	<p>Lower Extremity Angle Ranges from 0 to 78°</p> <p>Head Section Angle Ranges from +45° to -90°</p>	<p>Lower Extremity Angle Ranges from 0 to 78°</p> <p>Head Section Angle Ranges from +45° to -90°</p>	<p>Head Section Angle Ranges from +45° to -90°</p> <p>Lower Extremity Angle Ranges from 0 to 78°</p>
AM-134	AM-234	AM-334	AM-368	AM-368PB

OPTIONS:

<input type="radio"/> AM-1028 - Add Base with 4 Total Lock Swivel Casters (34"W Only)
<input type="radio"/> AM-1064 - Add Base with 4 Casters and 6" Clearance (34"W Only, Adds 4")
<input type="radio"/> AM-17 - Face Plug
<input type="radio"/> AM-91 - Add 1" Extra Soft Foam
<input type="radio"/> AM-807 - Rail Side Fold Down (Pair) - 34"W Only
<input type="radio"/> 19859 - Add-on Belt & Brackets
<input type="radio"/> SPCO - Special Vinyl Color (Provided by Customer)
<input type="radio"/> SPBASE - Special Base Color (Select Custom Base Color)
<input type="radio"/> 12088 - Side Panel 9 x 33.25 Perforated
<input type="radio"/> 12075 - Side Panel 9 x 33.25 Flat




STANDARD UPHOLSTERY CHOICES:

<input type="radio"/>  Imperial Blue	<input type="radio"/>  Taupe
<input type="radio"/>  Blue Ridge	<input type="radio"/>  Burgundy
<input type="radio"/>  Forest Green	<input type="radio"/>  Dove Gray
<input type="radio"/>  Black	<input type="radio"/>  Greystone
<input type="radio"/>  Cappuccino	<input type="radio"/>  Merlot
<input type="radio"/>  Patina	

STANDARD BASE COLOR:

<input type="radio"/>  Office White
--

CUSTOM BASE COLORS:

<input type="radio"/>  Black
<input type="radio"/>  Steel Gray
<input type="radio"/>  Metallic Silver