

■ SELECT COLLECTION

# Siesta Key Sleeper<sup>®</sup>

Available: Manual

Safety. Comfort. Cleanliness. Ambiance. Every aspect of the labor, delivery, and birthing experience can make a lasting impression on new mothers and their families.

The open-concept design of the Siesta Key Sleeper provides stylish daytime seating and spacious sleeping for overnight guests.

This sleeper features time-tested engineering with superior durability for years of use.



## Key Benefits

- Three functional configurations make product perfect for patient rooms and overnight visitor areas.
- Easy operation, low maintenance, comfortable. Powder-coated steel frame construction for durability, and rear-locking caster wheels for mobility.
- For quality engineering and patient-centric innovation that prioritizes functionality, safety, and comfort, Champion products comply with ISO 13485 and are made in the USA.

## Common Applications:

Labor & Delivery Rooms, General Patient Room Seating

## Standard Features



Steel and Wood Construction



Multiple Positions




Open Design




Caster Wheels (Rear-Locking)

## Space Requirements


Front View




Top View



Side View



Top View laying flat



## Specifications

<b>Model Number</b>	526
<b>Weight Capacity</b>	300 lbs. / 136 kg
<b>Overall Height</b>	42" / 106.70 cm
<b>Back Height</b> (above seat)	25" / 63.50 cm
<b>Overall Width</b>	32.50" / 82.60 cm
<b>Sleep Length</b>	74" / 187.96 cm
<b>Seat Height</b> ♦ <b>Seat Width</b>	18.75" ♦ 25" 47.63 cm ♦ 63.50 cm
<b>Seat Depth</b>	21" / 53.34 cm
<b>Arm Width</b>	2.50" / 6.35 cm
<b>Total Number of Positions</b>	3