

AM-BAX SERIES

HI-LO TREATMENT TABLES

STANDARD FEATURES:

- Load capacity: 500 lbs
- 27"W x 72 or 76"L x 19-38"H
- 34 oz. Heavy-duty vinyl with Permablok 3® bacterial protection
- 1" High density foam
- Rounded corners
- (4) Total locking casters
- Steel base with dual foot bar control
- White powder coat base



AM-BAX1000

TABLE CONFIGURATIONS:

<input type="radio"/> AM-BAX1000 - 1 Section 27"W x 76"L x 19-38"H	<input type="radio"/> AM-BAX2000 - 2 Section 27"W x 72"L x 19-38"H	<input type="radio"/> AM-BAX2500 - 2 Section with 3pc Head Section 27"W x 72"L x 19-38"H	<input type="radio"/> AM-BAX3500 - 3 Section 27"W x 72"L x 19-38"H
QTY	QTY	QTY	QTY
<input type="radio"/> AM-BAX5400 - 3 Section with Crescent Armrest 27"W x 72"L x 19-38"H	<input type="radio"/> AM-BAX3000 - 3 Section with Electric Elevating Midsection 27"W x 72"L x 19-38"H	<input type="radio"/> AM-BAX5500 - 5 Section with Adjustable Armrests 27"W x 72"L x 19-38"H	<input type="radio"/> AM-BAX5000 - 5 Section with Electric Elevating Midsection & Adjustable Armrests 27"W x 72"L x 19-38"H
QTY	QTY	QTY	QTY

TABLE FEATURES:

 <p>Flat Cushion</p>	 <p>Head Section Angle Ranges from +25° to -80°</p>	 <p>Head Section Angle Ranges from +25° to -45°</p>	 <p>Head Section Angle Ranges from +25° to -80°</p>
AM-BAX1000	AM-BAX2000	AM-BAX2500	AM-BAX3500
 <p>Lower Extremity Angle Ranges from 0 to 75°</p> <p>Head Section Angle Ranges from +25° to -45°</p>	 <p>Electric Elevating Midsection</p> <p>Lower Extremity Angle Ranges from 0 to 73°</p> <p>Head Section Angle Ranges from +45° to -90°</p>	 <p>Lower Extremity Angle Ranges from 0 to 73°</p> <p>Head Section Angle Ranges from +45° to -90°</p> <p>Armrests Vertical Adjustment 8" down</p>	 <p>Armrests Vertical Adjustment 8" down</p> <p>Electric Elevating Midsection</p> <p>Head Section Angle Ranges from +45° to -90°</p>
AM-BAX5400	AM-BAX3000	AM-BAX5500	AM-BAX5000

OPTIONS:

<input type="radio"/> AM-27 - Face Plug
<input type="radio"/> AM-91 - Add 1" Extra Soft Foam
<input type="radio"/> 19855 - Add-on Belt & Brackets
<input type="radio"/> SPCO - Special Vinyl Color (Provided by Customer)
<input type="radio"/> SPBASE - Special Base Color (Select Custom Base Color)

STANDARD BASE COLOR:

<input type="radio"/>	Office White
-----------------------	--------------

CUSTOM BASE COLORS:

<input type="radio"/>	Black
<input type="radio"/>	Steel Gray
<input type="radio"/>	Metallic Silver

STANDARD UPHOLSTERY CHOICES:

<input type="radio"/>	Imperial Blue	<input type="radio"/>	Taupe
<input type="radio"/>	Blue Ridge	<input type="radio"/>	Burgundy
<input type="radio"/>	Forest Green	<input type="radio"/>	Dove Gray
<input type="radio"/>	Black	<input type="radio"/>	Greystone
<input type="radio"/>	Cappuccino	<input type="radio"/>	Merlot
<input type="radio"/>	Patina		