

AM-SP SERIES

HI-LO TREATMENT TABLES

STANDARD FEATURES:

- Load capacity: 400 lbs
- 27"W x 76"L x 24-35.5"H
- One-step caster system
- Single pedestal base
- 34 oz. Heavy-duty vinyl with Permablok 3® bacterial protection
- 1.5" High density foam
- Rounded corners
- Single foot control pedal
- White powder coat base



NOTE: International (240V) option available

TABLE CONFIGURATIONS:









<input type="radio"/> AM-SP100 - 1 Section	<input type="radio"/> AM-SP200 - 2 Section	<input type="radio"/> AM-SP202 - 2 Section with Prenatal Cutout
QTY	QTY	QTY
		
<input type="radio"/> AM-SP350 - 3 Section	<input type="radio"/> AM-SP300 - 3 Section with Electric Elevating Midsection	<input type="radio"/> AM-SP500 - 5 Section with Electric Elevating Midsection
QTY	QTY	QTY
		

TABLE FEATURES:

 <p>Flat Cushion</p>	 <p>Head Section Angle Ranges from +40° to -90°</p>	 <p>Lower Extremity Angle Ranges from 0 to 75°</p>
AM-SP100	AM-SP200	AM-SP202
 <p>Lower Extremity Angle Ranges from 0 to 73°</p> <p>Head Section Angle Ranges from +45° to -90°</p>	 <p>Lower Extremity Angle Ranges from 0 to 73°</p> <p>Head Section Angle Ranges from +45° to -90°</p> <p>Electric Elevating Midsection</p>	 <p>Lower Extremity Angle Ranges from 0 to 73°</p> <p>Head Section Angle Ranges from +45° to -90°</p> <p>Armrests Range from +45° to -90°</p>
AM-SP350	AM-SP300	AM-SP500

OPTIONS:

<input type="radio"/> AM-SP1024 - Add Base with 4 Total Lock Swivel Casters
<input type="radio"/> AM-18 - Face Cutout
<input type="radio"/> AM-17 - Face Plug (Must order AM-18 also if ordering AM-SP100)
<input type="radio"/> AM-91 - Add 1" Extra Soft Foam
<input type="radio"/> SP30 - 30"W
<input type="radio"/> 19859 - Add-on Belt & Brackets
<input type="radio"/> SPCO - Special Vinyl Color (Provided by Customer)
<input type="radio"/> SPBASE - Special Base Color (Select Custom Base Color)

STANDARD BASE COLOR:

<input type="radio"/>	Office White
-----------------------	--------------

CUSTOM BASE COLORS:

<input type="radio"/>	Black
<input type="radio"/>	Steel Gray
<input type="radio"/>	Metallic Silver

STANDARD UPHOLSTERY CHOICES:

<input type="radio"/>	Imperial Blue	<input type="radio"/>	Taupe
<input type="radio"/>	Blue Ridge	<input type="radio"/>	Burgundy
<input type="radio"/>	Forest Green	<input type="radio"/>	Dove Gray
<input type="radio"/>	Black	<input type="radio"/>	Greystone
<input type="radio"/>	Cappuccino	<input type="radio"/>	Merlot
<input type="radio"/>	Patina		