

AM-SX SERIES

HI-LO TREATMENT TABLES

STANDARD FEATURES:

- Load capacity: 500 lbs
- 27"W x 72"L x 18-37"H
- 34 oz. Heavy-duty vinyl with Permablok 3® bacterial protection
- 1" High density foam
- Rounded corners
- (4) Total locking casters
- White powder coat base



AM-SX1000

TABLE CONFIGURATIONS:






<input type="radio"/> AM-SX1000 - 1 Section	<input type="radio"/> AM-SX2000 - 2 Section	<input type="radio"/> AM-SX3500 - 3 Section
QTY	QTY	QTY
		
<input type="radio"/> AM-SX3000 - 3 Section with Electric Elevating Midsection	<input type="radio"/> AM-SX5000 - 3 Section with Electric Elevating Midsection & Adjustable Armrests	
QTY	QTY	
		

TABLE FEATURES:

 <p>Flat Cushion</p>	 <p>Head Section Angle Ranges from +25° to -80°</p>	 <p>Lower Extremity Angle Ranges from 0 to 78°</p> <p>Head Section Angle Ranges from +25° to -80°</p>
<p>AM-SX1000</p>	<p>AM-SX2000</p>	<p>AM-SX3500</p>
 <p>Lower Extremity Angle Ranges from 0 to 78°</p> <p>Head Section Angle Ranges from +25° to -80°</p> <p>Electric Elevating Midsection</p>	 <p>Lower Extremity Angle Ranges from 0 to 78°</p> <p>Head Section Angle Ranges from +25° to -80°</p> <p>Electric Elevating Midsection</p> <p>Armrests Vertical Adjustment 8" down</p>	
<p>AM-SX3000</p>	<p>AM-SX5000</p>	

OPTIONS:

- AM-SX60** - Add Base with Casters & 6" Clearance
- AM-27** - Face Plug
- AM-91** - Add 1" Extra Soft Foam
- 19855** - Add-on Belt & Brackets
- SPCO** - Special Vinyl Color (Provided by Customer)
- SPBASE** - Special Base Color (Select Custom Base Color)

STANDARD BASE COLOR:

-

CUSTOM BASE COLORS:

-
-
-

STANDARD UPHOLSTERY CHOICES:

- | | |
|-------------------------------------|---------------------------------|
| <input type="radio"/> Imperial Blue | <input type="radio"/> Taupe |
| <input type="radio"/> Blue Ridge | <input type="radio"/> Burgundy |
| <input type="radio"/> Forest Green | <input type="radio"/> Dove Gray |
| <input type="radio"/> Black | <input type="radio"/> Greystone |
| <input type="radio"/> Cappuccino | <input type="radio"/> Merlot |
| <input type="radio"/> Patina | |