


HOW TO MEASURE - FLEXIBACK FRONTAL



GENERAL MEASUREMENTS

1. CHEST

Measure the circumference under your arms around the widest part of your chest.

2. WAIST

Measure the circumference around the narrowest part of your waist.

3. HIPS

Measure the circumference around the widest part of your hips, keeping your feet close together.

ONE PIECE APRON

Measure the length from the highest point of your trapezius to the top of your knee.

SIZING CHART

UNISEX

| Size Name | Chest Size | Length |
|-----------|------------|--------|
| XS | 32" - 34" | 30" |
| Small | 34" - 36" | 32" |
| Medium | 36" - 38" | 34" |
| Large | 38" - 40" | 36" |
| XL | 40" - 43" | 38" |
| 2XL | 43" - 46" | 40" |
| 3XL | 46" - 50" | 42" |