


TruAir Scope Drying Cabinets

Model Numbers:	SAS-18P-DRY, SAS-9P-DRY, SAS-8P-DRY, SAS-4P-DRY	Ext. Dimensions:	18-3/4"D X 40-7/8"W X 93-3/4"H (SAS-18P-DRY) 18-3/4"D X 40-7/8"W X 93-3/4"H (SAS-9P-DRY) 18-3/4"D X 21-5/8"W X 93-3/4"H (SAS-8P-DRY) 18-3/4"D X 21-5/8"W X 93-3/4"H (SAS-4P-DRY)
Color:	Bright White	Construction:	Powder Coated Metal Locking Hinged Glass Doors
Air Supply:	Air Pump (Upper Manifold) HEPA Filter Fan (Cabinet)	Air Flow Rate Filter Efficiency:	6 L/min per Scope (Minimum)
Air Pressure:	4 psi (Maximum)	Noise:	45 dB
Linear Pump Filter Efficiency:	99.97% up to 0.3µm	Cabinet Filter Efficiency:	99.97% up to 0.3µm
Electrical Supply:	120 V Wall Plug-3 A (10' Cord Exiting Top Right of Cabinet)	Power Consumption:	128W
Rated Amperage:	4.1A	Capacity:	18 Flexible Endoscopes


Measurements



POWDER-COATED METAL CABINET
 SAS-18P-DRY Capsa METAL SCOPE CABINET
 Locking Overlay Hinged Doors with Glass Inserts & Sloped Top, Holds 18 Flexible Endoscopes

TruAir Scope Drying Cabinets


Measurements



POWDER-COATED METAL CABINET
 SAS-9P-DRY Capsa METAL SCOPE CABINET
 Locking Overlay Hinged Doors with Glass Inserts & Sloped Top, Holds 9 Flexible Endoscopes

TruAir Scope Drying Cabinets


Measurements



POWDER-COATED METAL CABINET
SAS-8P-DRY Capsa METAL SCOPE CABINET
 Locking Overlay Hinged Doors with Glass Inserts & Sloped Top, Holds 8 Flexible Endoscopes

TruAir Scope Drying Cabinets

Measurements



POWDER-COATED METAL CABINET
 SAS-4P-DRY Capsa METAL SCOPE CABINET
 Locking Overlay Hinged Doors with Glass Inserts & Sloped Top, Holds 4 Flexible Endoscopes