




CAPSAHEALTHCARE



Secure-A-Scope™ Endoscope Storage Cabinet


Model Numbers:	SAS-4004, SAS-4008 (Half Size) SAS-4009, SAS-4018 (Full Size)	Ext. Dimensions:	93-3/4"H X 21-5/8"W X 18-3/4"D (Half Size) 93-3/4"H X 40-7/8"W X 18-3/4"D (Full Size)
Color:	Bright White	Construction:	Powder Coated Metal Locking Hinged Glass Doors Padded Lower Walls (Distal Tip Protection)
Cabinet Filter Efficiency:	Removes 99.97% of Particles (up to 0.3µm)	Air Supply:	HEPA Filter Fan (Cabinet)
Capacity:	4, 8, 9, or 18 Flexible Large Diameter Endoscopes	Electrical Supply:	120 V Wall Plug-3 A (8' Cord Exiting Top Right of Cabinet)
		Options:	Door Activated LED Lighting Keyless Lock Temperature and Humidity Gauge Stainless Steel

Measurements



Item
POWDER-COATED METAL CABINET
SAS-4004 Capsa METAL SCOPE CABINET
Locking Overlay Hinged Door with Glass Insert, Holds 4 Flexible Endoscopes

Secure-A-Scope™ Endoscope Storage Cabinet



Item
POWDER-COATED METAL CABINET
SAS-4008 Capsa METAL SCOPE CABINET
 Locking Overlay Hinged Door with Glass Insert, Holds 8 Flexible Endoscopes




CAPSAHEALTHCARE



Secure-A-Scope™ Endoscope Storage Cabinet

Measurements



Item
 POWDER-COATED METAL CABINET
 SAS-4009 Capsa METAL SCOPE CABINET
 Locking Overlay Hinged Doors with Glass Inserts & Sloped Top, Holds 9 Flexible Endoscopes




CAPSAHEALTHCARE



Secure-A-Scope™ Endoscope Storage Cabinet

Measurements



Item
 POWDER-COATED METAL CABINET
 SAS-4018 Capsa METAL SCOPE CABINET
 Locking Overlay Hinged Doors with Glass Inserts & Sloped Top, Holds 18 Flexible Endoscopes