



1	FLAT TOP W/ INVERTED "V" EDGE 4X
--	16GA300SERIES SS
2	1" DIAM. SS TUBE HANDLE
3	1" SS FLAT BAR HANDLE BRKT.
4	1-1/2" SQ. SS TUBE LEGS
5	SLIDES FOR BINS W/RAISED ENTRY
--	16GA300SERIES SS
6	HIGH TEMP PHENOLIC CASTERS 4" Dia.
7	ENDOSCOPE BINS



Item
 STAINLESS STEEL CART
 M6204SS MASS STAINLESS STEEL, STEEL CLEANABLE
 (Autoclavable) 4-Tier Scope Transport Cart