

accucold®

Medical-Life Science



DL11BWIFIKIT



**Wi-Fi & Bluetooth Enabled
For Real-Time Monitoring**



**Stores up to
32,000 Readings**



**Audible Alarm & Red LED
Excursion Signals**



**ISO/IEC 17025 NIST Calibration
Certificate Included**



**User-Friendly By Design
Free Intuitive Software**



Smart Probe Ready



**Secure Data Collection
Meets: FDA 21 CFR Part 11**

Cloud Monitoring Solution provides instant data access and alert notifications to ensure that you never miss a critical event.

accucold®

Memory Capacity	32,000 Total Memory 16,000 Per Channel - (Optional 2nd Channel)
Subscription	Free Tier Available 2yr Premium Account - \$95 - Email, Text, or Call Options!
Battery	3 Power Sources 5V Micro USB (Main); 2 x AAA Battery; & 1 x CR2032 Battery (backup)
Buffered Probe	Non-Breakable Bottle with Glass Beads
Alerts	Receive notifications as SMS, Email, WhatsApp, or Push notifications on your device.
Real-Time Data	Data is uploaded to the cloud in real time.
Accessibility	Download data files and reports onto your device.


The **Mobile App** enables users to access their **Premium Cloud** Account via an application, allowing them to record, store, and share logger data, receive alarm notifications, and view logger information without requiring a computer.



Calibrated Smart Probe



Buffer Bottle



Wall Mount



Wall Plug



Micro USB Cable



AAA Batteries

Single Channel SKU : DL11BWIFIKIT
Dual Channel SKU : DL22BWIFIKIT