

LockBoxWKit



5.88" x 19.13" x 10.25" (H x W x D)

Interior locking compartment with combination lock and front windows for use with select Accucold and Summit refrigerators

Highlights:

Designed for installation inside your refrigerator for an added level of content security

Manufactured in the USA from powder coated steel

Two windows provide a clear view inside

Product Features:

Combination lock	User-set combination lock allows you to change the code as needed to best accommodate your security needs
Key override	Two keys included to override the combination if needed
Powder coated steel constructed	Designed and manufactured in the USA
Windows included	Includes two clear plexiglass windows to provide a convenient view of stored contents; available without windows by special order
Flexible design	Works with a variety of Accucold and Summit refrigerators

