

TM1AKIT



4.25" x 5" x 1.75" (H x W x D)

Remote alarm contacts for aftermarket use on most refrigeration equipment with a display of the current and high/low temperature and audible/visual temperature alarm

Highlights:

Dry contacts with remote alarms for easy connection to your facility's monitoring system

Large LCD screen keeps the current and min/max temperatures in view at all times

High/low temperature alarm provides audio and visual warning if temperature rises or falls out of user-set range

Product Features:

Remote alarm contacts	Includes contacts to forward alarms to external remote system
Compatible on any Accucold or Summit model	Can be installed on virtually any brand of refrigerator or freezer. To have this factory-installed on an Accucold or Summit unit, add upgrade "TM1" when ordering.
Temperature alarm	Audible alarm will sound and indicator light will flash if temperature rises or falls below the specified approved range
Alarm delay	Allows you to delay alarm activation by a user-set time interval to avoid unnecessary alerts
Mute function	Selectable alarm mute period to best meet your needs
+/-1°C accuracy	Reliable temperature monitoring when used within a temperature range of -100 to +100°C
Protected temperature sensor	Temperature sensor is enclosed inside glycol-filled bottle to protect it from being affected by rapid temperature changes that may occur when the unit's door is opened
Extended sensor cable	10' long Teflon sensor cable for the temperature probe allows you to set the probe in your preferred location
Large simultaneous temperature display	LCD screen keeps the current, minimum, and maximum temperatures in view at all times, with a readout to the nearest tenth of a degree

TM1AKIT Specifications:

Overview	
Nominal Height	4.25" (11 cm)
Height of Cabinet	4.25" (11 cm)
Nominal Width	5.0" (13 cm)
Width	5.0" (13 cm)
Nominal Depth	2.0" (5 cm)
Depth	1.75" (4 cm)
Cabinet	Blue
Other Energy Info	12VDC Power Supply
Weight	1.0 lbs. (0 kg)
Parts & Labor Warranty	1 Year
UPC	761101103501

