

# ACR1151



58.38" x 23.63" x 22.88" (H x W x D)

Mid-sized pharmaceutical all-refrigerator with stainless steel construction inside and out, digital controls, and self-closing glass door

### Highlights:

Stainless steel construction inside and out for easier cleaning and improved sanitation

Self-closing door and open door alarm ensure items are not accidentally overexposed to warm temperatures

Full display all the way down to the floor

### Product Features:

|                                              |                                                                                                                        |
|----------------------------------------------|------------------------------------------------------------------------------------------------------------------------|
| Stainless steel interior and exterior finish | Cabinet and interior walls are constructed from stainless steel                                                        |
| Frost-free                                   | No-frost operation ensures minimum user maintenance                                                                    |
| Self-closing door                            | Door will automatically close when left ajar, preventing items from accidentally being exposed to ambient temperatures |
| Glass door                                   | Tempered glass door with stainless steel trim provides heat-safe view of interior contents                             |
| Recessed LED lighting                        | Energy efficient lighting offers full view of interior                                                                 |
| Factory installed lock                       | Keyed lock for a secure interior                                                                                       |
| Commercial design                            | Meets UL-471 commercial standards                                                                                      |
| Adjustable cantilevered shelves              | Adjustable shelving accommodates virtually any size with convenient cantilevered system                                |
| Fan-cooled interior                          | Internal fan ensures the interior is evenly cooled from top to bottom                                                  |

**ACR1151 Specifications:**

| Overview                   |                       |
|----------------------------|-----------------------|
| Nominal Height             | 59.0" (150 cm)        |
| Height of Cabinet          | 58.38" (148 cm)       |
| Nominal Width              | 24.0" (61 cm)         |
| Width                      | 23.63" (60 cm)        |
| Nominal Depth              | 23.0" (58 cm)         |
| Depth                      | 22.88" (58 cm)        |
| Depth with Handle          | 24.0" (61 cm)         |
| Depth with door at 90°     | 45.0" (114 cm)        |
| Capacity                   | 11.0 cu.ft. (311 L)   |
| Defrost Type               | Frost-Free            |
| Door                       | Glass                 |
| Cabinet                    | Stainless Steel       |
| US Electrical Safety       | ETL                   |
| Canadian Electrical Safety | ETL-C                 |
| Amps                       | 1.7                   |
| Voltage/Frequency          | 115 V AC/60 Hz        |
| Weight                     | 141.0 lbs. (64 kg)    |
| Parts & Labor Warranty     | 1 Year                |
| Compressor Warranty        | 5 Years               |
| UPC                        | 761101050348          |
| Refrigerator Features      |                       |
| Door Swing                 | RHD                   |
| Reversible                 | No                    |
| Shelf Type                 | Chrome                |
| Shelf Qty                  | 6                     |
| Adjustable Shelves         | Yes                   |
| Thermostat Type            | Digital               |
| Fan Type                   | Interior and Exterior |
| Refrigerant Type           | R600a                 |
| Refrigerant Amount         | 1.66oz.               |
| High Side PSI              | 270.0                 |
| Low Side PSI               | 105.0                 |
| Level Legs Qty             | 4                     |
| Interior Light             | Yes                   |
| Temperature Range          | 2° to 8°C             |
| Dimensions                 |                       |
| Interior Height            | 50.25" (128 cm)       |

|                        |                |
|------------------------|----------------|
| Interior Width         | 19.38" (49 cm) |
| Interior Depth         | 17.75" (45 cm) |
| Compressor Step Height | 4.75" (12 cm)  |
| Compressor Step Width  | 19.38" (49 cm) |
| Compressor Step Depth  | 5.75" (15 cm)  |

