

**CERVICAL PATIENT POSITIONING SANDBAG
7 LB 10" X 10"**

PRODUCT INFORMATION

PRODUCT INFORMATION

CERVICAL PATIENT POSITIONING SANDBAG - 7 LB 10" X 10"



The Cervical Patient Positioning Sandbag - 7 lb, 10" x 10" is crafted from premium, medical-grade heavy-duty black vinyl to provide dependable weighted immobilization for x-ray procedures. Designed to help technologists position patients with accuracy and ease, each sandbag features a heat-sealed inner sleeve to prevent leakage, a reinforced outer cover for long-lasting durability, and sturdy webbing handles for easy placement and transport. The smooth vinyl surface cleans quickly with any medical-grade disinfectant and supports use across therapeutic, rehabilitation, interventional, general positioning, and veterinary applications.

SPECIFICATION

COLOR: Black

MATERIAL: Vinyl

SIZE: 10" x 10"

WEIGHT: 7 lbs

FEATURES: Heat-sealed inner sleeve, Reinforced sewn and cover, Sturdy webbing handles, Easy to clean.

APPLICATIONS: Medical, Veterinary use, Therapeutic and rehabilitation purposes, Interventional after-care, General patient positioning.