

**CERVICAL PATIENT POSITIONING SANDBAG KIT**  
**5 LB, 7 LB, 10 LB**

**PRODUCT INFORMATION**

## PRODUCT INFORMATION

CERVICAL PATIENT POSITIONING SANDBAG KIT - 5 LB, 7 LB, 10 LB



The Cervical Patient Positioning Sandbag Kit - 5 lb, 7 lb, 10 lb is crafted from premium, medical-grade heavy-duty black vinyl to provide dependable weighted immobilization for x-ray procedures. Designed to help technologists position patients with accuracy and ease, each sandbag features a heat-sealed inner sleeve to prevent leakage, a reinforced outer cover for long-lasting durability, and sturdy webbing handles for easy placement and transport. The smooth vinyl surface cleans quickly with any medical-grade disinfectant and supports use across therapeutic, rehabilitation, interventional, general positioning, and veterinary applications.

This Cervical Patient Positioning Sandbag Kit includes six sandbags: two 5 lb bags (9" x 9"), two 7 lb bags (10" x 10"), and two 10 lb bags (11" x 11").

## SPECIFICATION

**COLOR:** Black

**MATERIAL:** Vinyl

**WEIGHT AND DIMENSIONS:**

- 2 x 5 lb. (9"x9")
- 2 x 7 lb. (10"x10")
- 2 x 10 lb. (11"x11")

**FEATURES:** Heat-sealed inner sleeve, Reinforced sewn and cover, Sturdy webbing handles, Easy to clean.

**APPLICATIONS:** Medical, Veterinary use, Therapeutic and rehabilitation purposes, Interventional after-care, General patient positioning.