

**CERVICAL PATIENT POSITIONING SANDBAG SET  
5 LB 9" X 9"**

**PRODUCT INFORMATION**

## PRODUCT INFORMATION

CERVICAL PATIENT POSITIONING SANDBAG SET - 5 LB 9" X 9"



The Cervical Patient Positioning Sandbag Set - 5 lb, 9" x 9" is crafted from premium, medical-grade heavy-duty black vinyl to provide dependable weighted immobilization for x-ray procedures. Designed to help technologists position patients with accuracy and ease, each sandbag features a heat-sealed inner sleeve to prevent leakage, a reinforced outer cover for long-lasting durability, and sturdy webbing handles for easy placement and transport. The smooth vinyl surface cleans quickly with any medical-grade disinfectant and supports use across therapeutic, rehabilitation, interventional, general positioning, and veterinary applications.

This Cervical Patient Positioning Sandbag Set includes two sandbags, each weighing 5 lb and measuring 9" x 9".

## SPECIFICATION

---

**COLOR:** Black

**MATERIAL:** Vinyl

**SIZE:** 9" x 9"

**WEIGHT:** 5 lbs each

**FEATURES:** Heat-sealed inner sleeve, Reinforced sewn and cover, Sturdy webbing handles, Easy to clean.

**APPLICATIONS:** Medical, Veterinary use, Therapeutic and rehabilitation purposes, Interventional after-care, General patient positioning.